

---

# Use of Controlled-Release Silver Dressing (CRSD)\* for a Full-Thickness Wound on Irradiated Tissue

Maeve Curran, PT and Oscar J. Paz-Altschul, MD  
Desert Regional Medical Center, Palm Springs, California

## Introduction

Radiation therapy (RT) is used to treat a variety of diseases. Radiation works by damaging the DNA of the cell, inhibiting its ability to replicate. The patient is treated with a radiation beam that can leave him with full-thickness tissue damage including, dry, frail, hairless skin with an impaired ability to heal. This can result in disfigurement and dysfunction, especially on an extremity. RT is nonselective, meaning that any cell and tissue in the radiation field may be affected. Often, normal function of fibroblasts, white and red blood cells, blood vessels, sebaceous glands, sweat glands and hair follicles will be affected. Patients should be taught to protect this area from further injury and trauma as they are at risk for developing wound site infections. The new advances in wound care protect the patient from pathogenic organisms that can lead to infection, delayed healing and provide the optimal moist environment needed for healing.

\* *Arglaes from Medline Industries, Inc.*

## **Material and Methods**

This is a 60 year old male that presented with a squamous cell carcinoma on the dorsum of his left hand. He underwent RT for five weeks resulting in a full-thickness wound. He was in his usual state of good health and nutrition was adequate. His lab values were within normal limits.

## **Discussion**

The proven antimicrobial effects of silver are well documented in the literature. Historically, silver nitrate and silver sulfadiazine have enjoyed popularity in the treatment of burns and wound care. There is no known resistance to silver and the CRSD\* has the ability to deliver ionic silver constantly and consistently for up to seven days, overcoming some of the negative aspects of these prior silver delivery methods. This CRSD has proven to be effective against a broad range of microorganisms including resistant strains such as MRSA and VRE. We chose this dressing for our patient because of its ability to deliver a controlled amount of ionic silver while maintaining an optimal moist wound healing environment.

## **Results**

The wound was completely reepithelialized in 16 days, utilizing appropriate wound care. This included aggressive cleansing, irrigation, debridement (both sharp and enzymatic) and antimicrobial silver dressings.

## **Conclusions**

This case study demonstrates how proper assessment of the wound, and early aggressive cleansing, and debridement allowed this wound to progress through the normal stages of wound healing from inflammatory to proliferative and currently, in the maturation phase. This also demonstrates that by using silver technology with its known antimicrobial effect, this patient had an excellent outcome and complete resolution without any complications.


---

## Case Studies



### **Day 1 (12/12/01) Pre-Arglaes**

*Full thickness wound, measuring 6.5 cm x 8.2 cm x 0.1 cm, with 100% necrotic tissue. Wound culture was positive for staphylococcus. The wound was irrigated with normal saline, sharply debrided, a enzymatic ointment applied, covered with a sterile 4 x 4 gauze dressing, and secured.*



### **Day 2 (12/13/01) Initiation of Arglaes**

*Dressing change, wound irrigation, and further sharp debridement was scheduled, however, the wound had completely debrided, leaving a 100% granular, full thickness wound. At this time, due to the potential for infection, the controlled release silver dressing was applied.*



### **Day 9 (12/21/01)**

*The wound base remained with 100% granulation tissue, measuring 3 cm x 4 cm with no depth.*



### **Day 16 (12/28/01)**

*The wound was completely reepithelialized.*

## References

Bates-Jensen B. (2001) "Management of Necrotic Tissue." In Sussman, C., Bates-Jensen, B. (Eds.). *Wound Care: A Collaborative Practice Manual for Physical Therapists and Nurses, Second Edition*. Aspen Publishers, Inc., Gaithersburg, Maryland.

Dow G, Browne A and Sibbald R.G. (1999) "Infection in Chronic Wounds: Controversies in Diagnosis and Treatment." *Ostomy Wound Management*; 45(8): 23-7, 29-40.

Falconio-West M and Sacramento L. (June 1999) "Silver Dressing Improves Healing in MRSA Wound," *Presented at the 31st Annual Wound, Ostomy and Continence Conference*, Minneapolis, Minnesota.

Martinez J. (April 1999) "Use of Silver Impregnated Transparent Film Dressing as the Primary Treatment of Pyoderma Gangrenosum," *Presented at the Symposium on Advanced Wound Care and Medical Research Forum on Wound Repair*, Anaheim, California

Rudolph R. (1982) "Complications of Surgery for Radiotherapy," *Skin Damage. Plastic Reconstruction Surgery* 70:179.

