QUICK START GUIDE

HEART RATE WATCH

Item # MDSP3043

Getting Started

- The watch requires 1 CR2032 type battery (included).
- Set the date and time before using this watch (see below).

See the enclosed User Manual for further details and instructions.

LCD DISPLAY

Modes

This watch has 4 main modes and several sub-modes.

- Time Display Mode (Main Screen)
 - o 2 Alarms (AL1 and AL2)
 - o Stopwatch
 - o Countdown Timer
- Pedometer Mode
 - o Total Steps/Current Time
 - o Total Steps/Exercise Time
 - o Calories/Distance
- Step Memory Storage Mode
- Heart Rate Mode

Setting the Time/Date

1. While in Time Display Mode, hold down MODE button for about 3 seconds to enter Time/Date setting mode.
2. Press the UP and SET buttons to adjust the settings up and down. Hold either button down for two seconds to make rapid adjustments.
3. Press MODE to switch between setting items. Order is: 12/24 Hour Format > Hour > Minute > Second > Year > Month > Day > SIG (hourly chime on/off) > Exit

Using the Pedometer Mode

1. From Time Display Mode, press MODE button once to switch to Pedometer Mode.
2. Press UP button to switch between sub-modes.
3. In the Steps/Time Mode or Steps/Exercise Time Modes, press and hold MODE for 2 seconds to enter User Information Setup.
4. Press MODE to switch between setting items. Use UP and SET buttons to edit each individual item. Order is Metric/Imperial System > Weight > Distance between Steps > G-Sensor Sensitivity > Exit
5. Press and hold UP button for 3 seconds to lock in this mode and prevent accidental mode switch. The watch will now accurately track your steps and calorie-loss, as long as you are in Pedometer Mode.

Using the Heart Rate Mode

1. From Time Display Mode, press MODE button three times to switch to Heart Rate Mode.
 - a. NOTE: Screen will say "PULSE" at the top with a heart icon next to it.
2. Use both fingers to press down on the metal pulse sensors at the top and bottom of the watch face. Cover the sensors as much as possible, and exert pressure evenly.
3. The heart icon will flash when a heartbeat is detected - keep fingers firmly pressed to the sensors.
4. The heart rate value will display and update in real time as your fingers remain pressed to the sensors.

Please read the Manual for this watch prior to using. If you have any questions or are experiencing difficulties, please call the Medline Service Center at **1-800-MEDLINE** (1-800-633-5463). RH17POA

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