
# NO-TOUCH FOREHEAD THERMOMETER


## APPEARANCE & KEY FUNCTIONS

## INTRODUCTION


Thank you for choosing the infrared forehead thermometer. This innovative medical device relies on advanced infrared (IR) technology to measure temperature instantly. The thermometer measures the infrared heat generated by the surface of the skin over the vessel and its surrounding tissue. The thermometer then converts it into a temperature value shown on the LCD screen.

## LCD SCREEN


1. Temperature unit
2. Temperature display
3. Temperature scanning in progress
4. Forehead temperature indicator
5. Low battery indicator
6. Object surface temperature indicator
7. Smiling face indicator (Temperature falls within normal range)
8. High temperature indicator


## SELECT TEMPERATURE UNIT (°C or °F)


1. Slide the Mode Switch from 'off' to 'forehead' mode and then immediately press and hold the 'Scan' button.

2. Keep holding down the 'Scan' button for 3-4 seconds until you see the word 'Set' appear on the screen.


3. After the word 'Set' appears, press the 'Scan' button to toggle between C and F temperature units.


4. After you have selected the preferred temperature unit, press and hold 'Scan' for 3-4 seconds until the screen readies for a forehead measurement (you will see the forehead symbol appear). The next measurement will display in the selected unit.

## TIPS FOR PERFORMING A TEMPERATURE MEASUREMENT


1. To take a temperature, the thermometer must be held steady within 1 inch (2-3 cm) from the forehead. Sweat and other obstructions, such as hair, will affect the accuracy of the reading. Please make sure that the forehead is clear and dry before attempting a measurement.


2. When the measurement is complete, use the Mode Switch to turn the thermometer off.

## PERFORM A MEASUREMENT


### SURFACE MODE

**NOTE:** While in surface mode, the measurement will perform continuously if you press and hold the scan button without release.

### FOREHEAD MODE

## RESULT INDICATORS

Temperature falls within measurement range: 93.9°F to 100.2°F (34.4°C to 37.9°C) and the LCD backlight will display "Green".


Temperature falls within measurement range: 100.4°F to 108°F (38.0°C to 42.2°C) and the LCD backlight will display "Red".


Temperature falls outside the displayed range:  
· Hi-temperature  $\geq 108.1^\circ\text{F}$  (42.3°C)  
· Lo-temperature  $\leq 93.7^\circ\text{F}$  (34.3°C) is shown on the display with "Blue" LCD backlight.


## REPLACING THE BATTERY


1. Remove battery cover

2. Remove the old battery with battery hook

3. Insert a new battery

4. Reattach the battery cover

## TROUBLESHOOTING

Message	What it means	What to do
Err1	Environmental temperature is below or above system operation range.	Put the thermometer under operating temperature range, 50°F to 104°F (10°C to 40°C)
Err2		
Err4	Problem with the thermometer.	Review the instructions and re-start the measurement. If the steps do not work, contact Medline customer service.
	Battery is low and 'Lo' appears on LCD.	The meter is functional and the result remains accurate, but it is time to change the battery.
Err5	The battery can't provide enough power for a test. After 3 beeps, the meter will turn off automatically.	Replace the battery immediately.

## DETAILED INFORMATION

Body temperature can vary from one individual/person to the next. It also varies by location on the body and time of day. Below shows the statistical normal ranges from different temperature sites. Please keep in mind that temperatures measured from different sites should not be directly compared, even if they are taken at the same time.

Table \*1 Body Site Normal Temperature Range

Oral	1°F (0.6°C) or more above or below 98.6°F (37°C)
Rectal/Ear	0.5°F to 1°F (0.3°C to 0.6°C) higher than oral temperature
Axillary (armpit)	0.5°F to 1°F (0.3°C to 0.6°C) lower than oral temperature

\*1. Body Temperature at WebMD; website: <http://firstaid.webmd.com/body-temperature> retrieved at 2010 Jan 7.

## WARRANTY

Medline Industries, Inc. warrants this precision digital thermometer against any defects in material and workmanship for a period of one (1) year after the date of purchase. This guarantee excludes cases of freight damage, complete water immersion, tampering, clear abuse, misuse, or accidents. Medline Industries, Inc. will, at its discretion, repair or replace this thermometer during the warranty period without charge. Claims must be filed together with proof of purchase within this warranty period. Battery life and/or function is not included under this warranty. No representative or person is authorized to assume for us any liability in connection with the sale of the products of Medline Industries, Inc. This warranty gives you specific legal rights, and you may also have other rights that may vary from state to state. Prior to shipping any product, please contract either Medline Industries, Inc. at 1-800-MEDLINE or an authorized Medline sales representative with any warranty concerns.

## IMPORTANT SAFETY INSTRUCTIONS


1. Close supervision is necessary when using this device around children or persons needing special care.
2. Use the thermometer only for the intended use described in this manual.
3. Do not use the thermometer if it is not working properly, or if it is damaged.
4. Do not use accessories that are not supplied or recommended by the manufacturer.
5. Proper maintenance is essential to the longevity of your device. If you are concerned about the accuracy of the device, please contact the manufacturer.

## CAUTIONS & WARNINGS

1. As with any thermometer, proper technique is crucial to getting accurate temperature readings. **Please read these instructions thoroughly and carefully before using the thermometer.**
2. Always operate the thermometer in an environment with a temperature range of 50°F to 104°F (10°C to 40°C), and relative humidity less than 95%.
3. Always store the thermometer in a cool and dry place: temperatures between -4°F to 140°F (-20°C to 60°C); relative humidity less than 95%.
4. **Avoid direct sunlight and avoid dropping the thermometer.**
5. **Do not** disassemble the thermometer.
6. Basic safety precautions should always be observed.
7. This thermometer is not intended to be a substitution for a consultation with your physician.

## CARE & CLEANING

- The probe is not water-proof. Please wipe with a clean, dry cotton swab to clean the probe.
- The body of the thermometer is water-resistant. Use a soft and dry cloth to clean it. Do not use abrasive cleaners.
- Store the thermometer in a cool, dry location, free from dust and away from direct sunlight.

## SYMBOL INFORMATION

-  Consult instructions for use
-  IP22 Ingress protection rating
-  Manufacturer
-  SN Serial number
-  Caution, consult accompanying documents
-  Authorized representative in the European Community
-  Collection for electrical & electronic equipment
-  Type BF Equipment
-  Temperature limitation
-  Humidity limitation

## REFERENCE STANDARDS

Device Standard:

Device Corresponds to the requirements of the standard for infrared thermometers

ASTM E1965-98, EN ISO 80601-2-56:2012,

EN 60601-1-2:2007/AC:2010,

EN 60601-1:2006+A12:2014, EN 60601-1-6:2010

Electromagnetic Compatibility:

Device fulfills the stipulations of the standard EN 60601-1-2

The stipulations of EU-Directive 93/42/EEC for

Medical Device Class II devices have been fulfilled.

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Manufactured for/Fabriqu  pour/Fabricado para: Medline Industries, Inc., Three Lakes Drive, Northfield, IL 60093 USA.