MEDLINE

REF MDSTH1002

TALKING INFRARED THERMOMETER

Please read this instruction manual before use.

CONTENTS

INTRODUCTION	3
IMPORTANT CLEANING INFORMATION	3
CARE AND CLEANING	4
IMPORTANT SAFETY INFORMATION	4
BODY TEMPERATURE	6
PACKAGE COMPONENTS	6
PRESENTATION	7
BEFORE USE	9
HOW TO MEASURE FOREHEAD TEMPERATURE	10
HOW TO MEASURE EAR TEMPERATURE	10
HOW TO MEASURE OBJECT TEMPERATURE	12
TEMPERATURE TAKING HINTS	12
TEMPERATURE BACKLIGHT INDICATOR	13
CHECK THE MEMORY	13
VOICE PROMPT	14
DATE AND TIME SETTING	15
CELSIUS/FAHRENHEIT OPTION	15
REPLACING THE BATTERY	16
SPECIFICATIONS	16
CALIBRATION	18
TROUBLESHOOTING	18
NORMALIZED SYMBOLS	20
DISPOSAL	20
WARRANTY	21
IMPORTANT INFORMATION REGARDING ELECTROMAGNETIC COMPATIBILITY (EMC)	22


INTRODUCTION

This thermometer has been carefully developed for accurate, safe, and fast human body temperature measurements via the ear and forehead. It has also been similarly validated for object temperature measurements.

This thermometer uses an offset when taking forehead and ear measurements so that the reading reflects the body's core temperature.

When used in accordance to the user's manual, this thermometer can be easily maintained. Please note that the ear is an ideal site for taking a body temperature measurement. This is because body temperature is regulated by the hypothalamus, which shares the same blood supply as the tympanic (ear) membrane.

The thermometer is an ear canal IR thermometer. The forehead temperature is obtained based on the site offset between the ear and the forehead, and the forehead temperature is only used as a reference.


This thermometer is intended for home use by customer only. The measure site is ear canal.

INTENDED USE

The device is intended for the intermittent measurement and monitoring of human body temperature by consumers in the home. It is intended for use on people of all ages.

IMPORTANT CLEANING INFORMATION

- Earwax in the ear canal and lipid, oil, or sweat on the forehead may cause an inaccurate measurement. Please make sure that the ear canal and/or forehead is clean before attempting a measurement.
- If you clean the ear or forehead, please waiting for 5-10 minutes before measuring.
- Don't touch or blow the infrared sensor. A polluted or contaminated infrared sensor may cause inaccuracy. If you think you have contaminated the probe sensor, please clean it.
- To clean the sensor, gently wipe its surface with a cotton swab slightly moistened with alcohol and immediately wipe dry with a clean cotton swab. After cleaning, allow at least 20 minutes of drying time before taking temperatures.

CARE AND CLEANING

To ensure accurate measurements, it is very important to keep the probe tip clean and free of scratches. Finger prints, earwax or dirt will affect the accuracy of the thermometer. In order to get an accurate measurement, please clean the probe tip.

Gently wipe its surface with a cotton swab slightly moistened with alcohol and immediately wipe dry with a clean cotton swab.

Use a soft dry cloth to clean the body of the thermometer. Never clean the thermometer with an abrasive cleanser, thinner, benzene or submerge the thermometer into water or other liquids.

After cleaning, place the thermometer in the protective case.

Store it in a clean, dry place at room temperature. Never expose the thermometer to extreme temperatures, humidity, direct sunlight or shock.

IMPORTANT SAFETY INFORMATION

Proper technique is critical in order to obtain accurate temperatures. Therefore, in order to obtain an accurate measurement, please read this manual carefully before using this thermometer.

- Use of this thermometer is not intended to substitute professional consultation by your physician. Measurement results are for reference only. Contact your physician about any known or suspected health concerns.
- The Infrared sensor must be kept clean, dry, and undamaged at all times to ensure accurate measurements.
- Do not use the ear/forehead measurement mode for measurements other than for human body temperature.
- Avoid measuring the ear temperature if the ear is inflamed or infected.
- If the temperature of the storage area differs greatly from that of the measuring area, please allow time for the thermometer to equalize to the current environment's temperature. For great differences in environment temperatures, let the thermometer equalize for 30 minutes before use.

- Keep the thermometer out of the reach of children and pets to avoid inhalation or swallowing of small parts. Do not allow children to take their temperatures unattended. Children may not be able to use the device according to the instructions in this user manual. This thermometer is not a toy.
- Don't throw the thermometer battery into fire. The battery is likely to explode in fire.
- Don't store the thermometer in extreme temperatures or extreme humidity. Doing so may cause inaccuracy.
- Don't operate the thermometer under extreme temperatures or extreme humidity. Doing so may cause inaccurate measurement.
- Don't use this product in ear mode inside of a wet ear, such as after showering, swimming, or bathing. Doing so may damage the ear or thermometer.
- Don't use the thermometer if it is damaged in any way (for example, the infrared sensor is broken). The continuous use of a damaged unit may cause injury, improper results, or serious danger.
- Don't take apart, repair or change any parts of the unit at any time outside of replacing the battery.
- Don't use a mobile phone nearby when the thermometer is in use.
- If your thermometer will not be used on a regular basis, remove battery to prevent possible damage due to chemical leakage. If the battery does leak, remove carefully. Do not allow bare skin to touch the leaking fluid. • When using, do not touch the patient and the battery simultaneously.
- The materials (e.g., ABS) of expected contact with the patient have passed the ISO 10993-10 standards test, meaning that this device will not cause toxicity, allergies, or irritation from normal use. However, if you have an allergic reaction to any material, please stop the use of this thermometer immediately and consult your physician.
- Store the thermometer in a dry, clean place. Do not leave the product exposed to any chemical solvent, lint, dust, direct sunshine or high temperatures.
- If you have any problems with this device, such as setting up, maintaining or using the thermometer, please contact the customer service department of the manufacturer, Medline Industries, Inc.
- Please report any unexpected events or reactions while using this device to the manufacturer.
- Please avoid mechanical shock. Otherwise the thermometer will be damaged.

Caution: This thermometer should not be used adjacent to or stacked with other equipment. If said equipment must remain in close proximity to the thermometer being used, it should be monitored in order to verify normal operation during the time of measurement.

BODY TEMPERATURE

Normal body temperature has a range. A person's normal temperature range tends to decrease with age. The following table shows normal temperature ranges by age:

0 – 2 years	36.4 – 38.0 °C	97.5 – 100.4 °F
3 – 10 years	36.1 – 37.8 °C	97.0 – 100.0 °F
11 – 65 years	35.9 – 37.6 °C	96.6 – 99.7 °F
> 65 years	35.8 – 37.5 °C	96.4 – 99.5 °F


The range of what is normal varies from person to person and can be influenced by many factors, such as time of day, measuring from different sites, level of activity, medications, emotion, and so on.

It is recommended that you practice with the thermometer on yourself and family members when you are healthy. This way you know how the thermometer works, and will feel more confident with the measurements when the same person is ill.


PACKAGE COMPONENTS

DESCRIPTION	QUANTITY
• Quick Start Guide	1
• Manual	1
• Clinical IR thermometer	1
• 1.5V AAA Battery	2

PRESENTATION

LCD DISPLAY

BEFORE USE

1. **To achieve accurate measurements, it is very important to check that the infrared sensor is clean before each use.**

To clean the sensor, gently wipe its surface with a cotton swab slightly moistened with alcohol and immediately wipe dry with a clean cotton swab. After cleaning, allow at least 20 minutes of drying time before taking temperatures.

2. Please open the battery cover and move the insulating piece when you use the thermometer for the first time. After its removal, the thermometer will turn on automatically, and the LCD display will show similar to the bottom picture on the left-hand side of this page.

Notes on the LCD Screen:

- 1) In measuring mode and memory mode, the LCD display is normal. In waiting mode and setting mode, the LCD display is reversed.
- 2) In waiting mode, the LCD display circulates between time, date, and room temperature.

HOW TO MEASURE FOREHEAD TEMPERATURE


1. Put on the forehead adaptor before use. Place the thermometer in your hand as shown on the right, with your thumb on the  button.

Make sure the probe is clean before measuring.

2. Hold the thermometer so that the forehead adaptor is flush with the center of the forehead, midway between eyebrow and hairline. Once the device is touching the forehead as indicated, hold the thermometer steady and press  to start the measurement.

3. Keep the thermometer flat until you hear the temperature voice announce the completed reading (if the talking feature is turned on), or until you hear a beep. The back-light will also light up according to the temperature reading. When the measurement is complete, the thermometer may be removed from the head.

4. Read the temperature on the display.

HOW TO MEASURE EAR TEMPERATURE

1. Remove the forehead adaptor before use. Place the thermometer in your hand as shown with your thumb or forefinger on the  button.

Make sure the probe is clean before measuring.

2. Insert the probe carefully into the ear canal. Hold the thermometer so that the probe is going straight into the ear toward the eardrum. It is important that the probe sensor faces the eardrum.

Hold the outer ear and gently pull it toward the rear of the head to straighten the ear canal. Always insert the thermometer into the same ear with the same direction and depth.

3. Press the  to start the measurement. Keep the probe still until you hear the temperature voice announce the completed reading (if the talking feature is turned on), or until you hear a beep. The back-light will also light up according to the temperature reading.

4. When the measurement is complete, the thermometer may be removed from the ear. The LCD will flash for about 6 seconds. You can measure again after the ear sign stops flashing.

HOW TO MEASURE OBJECT TEMPERATURE

1. Remove the forehead adaptor before use. Under waiting mode or measuring mode, place the thermometer in your hand, as shown, with your fingers on  and .

Press  and  to enter TARGET temperature measuring mode. The LCD display will look like the illustration on the right with the word 'food'. Make sure you do not use this mode to measure body temperature.


2. Hold the thermometer so that the probe is positioned at a distance within 1cm from the liquid and/or surface you want to measure. Press  button to start the measurement.

3. The LCD will display the measured temperature. You can press the  and  button again at the same time or press [MEM/SET] to return to waiting mode. The thermometer will automatically return to waiting mode after 1 minute.

TEMPERATURE TAKING HINTS

1. A temperature measurement taken in the right ear may differ from a measurement taken in the left ear. Therefore, always take the temperature in the same ear.
2. External factors may influence ear and forehead temperatures, particularly when an individual has:
 - been lying on one ear or the other
 - had their ears covered
 - been exposed to very hot or very cold temperatures, or been recently swimming or bathing.In all of the above cases, wait 30 minutes prior to taking a measurement.
3. For persons wearing hearing aids or ear plugs, remove such devices and wait 30 minutes prior to taking an ear temperature.
4. If ear drops or other topical ear medications are being used, use the untreated ear for measurements. If both ears are receiving treatment, please use the forehead measurement method.

TEMPERATURE BACKLIGHT FEVER INDICATION

This thermometer has a fever alarm indicator light on the end of the thermometer opposite from the probe sensor. This light will indicate if the measurement detects a fever:

Green (Normal): 96.0°F - 99.3°F (35.6°C - 37.4°C)

Slow Blinking Red (Raised): 99.4°F - 100.3°F (37.5°C - 37.9°C)

Fast Blinking Red (High): 'FE' will also display on the LCD, as indicated in the picture at left. This is for results between 100.4°F - 109.3°F (38.0°C - 42.9°C).

HOW TO CHECK THE MEMORY

This thermometer can recall the past 9 measurements.

The current measurement is always stored in the last storage space. When all storage spaces have been occupied, the oldest measurement is automatically deleted from the memory.

Press [MEM/SET] to check the memory number, temperature and date. See figure on the left.

Press [MEM/SET] repeatedly to cycle through the temperatures in the memory.

Press and hold [MEM/SET] to check the memory data.


Press  and  at the same time to return the device to *Waiting Mode*.

Waiting mode
screen.

VOICE PROMPT

Press and hold **[TALK/OK]** for approximately 2 seconds to enable or disable the voice prompt function. The device will beep and display the selected voice option when programming is complete. To switch between languages: English, Spanish, French, German, Italian and Russian*, repeat the last step until you land on the preferred language. If turned on, press **[TALK/OK]** to prompt the device to tell you the room temperature when in waiting mode, and your measured temperature when using the ear or forehead method.

*Not all features are supported by all languages.

DATE AND TIME SETTING

Press and hold **[MEM/SET]** about 2 seconds to go into the setting mode.

Press **[MEM/SET]** to switch the hour mode (12 hr or 24 hr), hour, minute, year (2001–2099), month, day and LCD light (ON or OFF).

Press **[MEM/SET]** to increase the parameter you selected.

Press the **[TALK/OK]** to confirm the current selection, and move to the next settings option.

When the sleep mode is active, the thermometer will turn off automatically about 1 min or 3 min (which is factory default); when the sleep mode is deactivated, the thermometer remains on.

To save power, we recommend turning sleep mode on.

Press  and  at the same time to enter *waiting mode*.

CELSIUS/FAHRENHEIT OPTION

Press and hold **[TALK/OK]** button to switch between Celsius (°C) and Fahrenheit (°F) when the thermometer is in measuring or memory mode.

REPLACING THE BATTERY

When the low battery indicator  appears on the LCD, or the thermometer does not function at all, you should replace the batteries as soon as possible.


Replacing the battery:

1. Remove the battery cover.
2. Observing the correct polarity, insert the 2x AAA batteries.
3. Replace the battery cover.
4. Dispose of used batteries according to the manufacturer's instructions.

WARNING:

If the battery leaks and comes into contact with the skin or eyes, wash immediately with copious amounts of water. Batteries should be handled under adult supervision.

- Dispose of used batteries with care.
- Do not disassemble, recharge, or dispose of batteries in fire.
- Do not mix new and old batteries.
- Keep the batteries out of the reach of children.

SPECIFICATIONS

Model number:	MDSTH1002/TH1002D
Power supply:	DC3V (2xAAA)
Temperature sensing part:	IR Sensor
Body measuring range:	89.6°F - 109.3°F (32.0°C - 42.9°C)
	Laboratory accuracy:
	±0.4°F (±0.2°C) in range
	of 95.9°F - 107.6°F (35.5°C - 42.0°C)
	±0.5°F (±0.3°C) for other range
Object measuring range:	68.0°F - 140°F (20.0°C - 60.0°C)
Object measure accuracy:	±4.0°F (±2.0°C)
Room temperature range:	41.0°F - 139.8°F (5.0°C - 59.9°C)
Room temperature accuracy:	±2°F (±1°C)
Measure time:	1 - 2 Seconds
Measure interval:	9S (For Ear mode)
Reference body site:	core
Operating mode:	adjusted mode
Resolution:	0.1°F (0.1°C)

Forehead operating condition:	59.0°F - 104.0°F (15.0°C - 40.0°C)
Ear operating condition:	50.0°F - 104.0°F (10.0°C - 40.0°C)
Relative humidity:	15%RH - 85%RH
Atmosphere pressure:	700hPa - 1060hPa
Storage and transportation conditions:	-4°F - 131°F (-20°C - 55°C)
Relative humidity:	15%RH - 90%RH
Atmosphere pressure:	700hPa - 1060hPa
Outer dimension:	142.5mm (L) x 38.5mm (W) x 31mm(H)
Weight:	63.5g (without batteries)
Service life of the device:	2 years
Service life of the battery:	With a new battery, approximately 6 months based on the use frequency of 5 times/day.

Note: Design and specifications are subject to change without notice.



This infrared thermometer meets requirements established in ASTM Standard E 1965-98 for the thermometer system. Full responsibility for the conformance of the product to the standard is assumed by Shenzhen Dongdixin Technology Co., Ltd, Floor 1-2, No.3 Building, Fanshen Xusheng Industrial Estate Xilixiaobaimang 518108 Nanshan District, Shenzhen P. R. China.

ASTM laboratory accuracy requirements for the thermometer only in the display range of 98 °F to 102.2 °F (37 °C to 39 °C) for infrared thermometers is ± 0.4 °F (± 0.2 °C), whereas for mercury-in-glass and electronic thermometers, the requirement per ASTM Standards E 667-86 and E 1112-86 is ± 0.2 °F (± 0.1 °C).

CALIBRATION

The thermometer is initially calibrated at the time of manufacture. If the thermometer is used according to the instructions, periodic recalibration is not required. If at any time you question the accuracy of the measurement, please contact the retailer immediately. Don't attempt to modify or reassemble the thermometer.

TROUBLESHOOTING

Problem	Checklists	Countermeasures
No response / Automatic reset when pulling out the insulator	Battery dead? Incorrect polarity? Poor battery contact?	Change batteries Take out battery, reinsert battery correctly
	The measured temperature is lower than 89.6°F /32.0°C or higher than 109.3°F /42.9°C Please check the operation method	Follow steps in User Manual for proper measurement and probe sensor placement.
	Hardware problem	Contact the manufacturer or distributor
	Please check if the infrared sensor is clean	Clean the infrared sensor with cotton swab according to user manual

Check to make sure you are using the correct measurement method for the sensor placement

Ensure you have read the manual and know how to use the thermometer properly

The measurement does not appear to be accurate.

Wait 30 minutes in the same environment and try again after the patient and thermometer have stabilized

Please keep the thermometer and patient in the same environment for at least 30 minutes before using

Are you using the thermometer indoors?

Please take the measurement indoors

The measurement is not accurate or does not appear to be accurate

If you are holding the thermometer in your hand, you could be warming it up unnecessarily, which could affect the accuracy

Put the thermometer on a flat surface and let it cool down first

Please check if there is a low battery icon on LCD

Replace with new batteries

NORMALIZED SYMBOLS

Applied part of type BF (The applied part is the probe tip)

Disposal in accordance with Directive 2012/19/EU (WEEE)

Complies with the European Medical Device Directive (93/42/EEC) and amended by directive 2007/47/EC requirements. Notified body TUV Rheinland (0197).

The name and the address of the manufacturer

The name and the address of the Authorized EC-representative in Europe

Refer to Instruction Manual.

Manufacturing date is given by the SN number located inside the battery compartment. From left to right, the first digits represent area code and the next four digits is manufacture year and month.

Date of Manufacture

WARRANTY

One year warranty is available from purchasing date, excluding caused failures listed below:

1. Failure resulted in unauthorized disassemble and modification.
2. Failure resulted from an unexpected drop.
3. Failure resulted from operation in a way other than what is specified in the User's Manual.

DISPOSAL

Exhausted batteries do not belong in household waste. Dispose of batteries according to the current federal, state and local regulation. As a consumer, you are obligated by law to return the spent batteries.

IMPORTANT INFORMATION REGARDING ELECTROMAGNETIC COMPATIBILITY(EMC)

Important information regarding Electro Magnetic Compatibility (EMC)

With the increased number of electronic devices such as PC's and mobile (cellular) telephones, medical devices in use may be susceptible to electromagnetic interference from other devices. Electromagnetic interference may result in incorrect operation of the medical device and create a potentially unsafe situation. Medical devices should also not interfere with other devices.

In order to regulate the requirements for EMC (Electro Magnetic Compatibility) with the aim to prevent unsafe product situations, the EN60601-1-2 standard has been implemented. This standard defines the levels of immunity to electromagnetic interferences as well as maximum levels of electromagnetic emissions for medical devices. MDSTH1002 conforms to this EN60601-1-2:2007 standard for both immunity and emissions.

Nevertheless, special precautions need to be observed:

- The medical devices should not be used adjacent to or stacked with other equipment. In case adjacent or stacked use is necessary, the medical device should be observed to verify normal normal operation in the configuration in which it will be used.
- Refer to further guidance below regarding the EMC environment in which the device should be used.

(Table 1)

Guidance and manufacturer's declaration —electromagnetic emissions		
IR Thermometers are intended for use in the electromagnetic environment specified below. The customer or users of this Clinical IR Thermometers should assure that it is used in such environments.		
Emissions test	Compliance	Electromagnetic environment guidance
RF emissions CISPR 11	Group 1	MDSTH1002 IR Thermometer uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	MDSTH1002 IR Thermometer is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

(Table 2)

Guidance and manufacturer's declaration –electromagnetic immunity			
MDSTH1002 IR Thermometer is intended for use in the electromagnetic environment specified below. The customer or users of this Clinical IR Thermometers should assure that it is used in such environments.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ±8 kV air	± 6 kV contact ±8 kV air	Floor should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC61000-4-4	Not applicable	Not applicable	Not applicable
Surge IEC 61000-4-5	Not applicable	Not applicable	Not applicable
Voltage dips short interruptions and voltage variations on power supply IEC 61000-4-11	Not applicable	Not applicable	Not applicable
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

(Table 4)

Guidance and manufacturer's declaration –electromagnetic immunity			
MDSTH1002 IR Thermometer is intended for use in the electromagnetic environment specified below. The customer or users of this Clinical IR Thermometers should assure that it is used in such environments.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	Not applicable	3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the MDSTH1002 IR Thermometer, including cables, than the recommended separation distance calculated from the equation appropriate to the frequency of the transmitter. Recommended separation distance $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 1.2 \sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters as determined by an electromagnetic site survey.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz		

(Table 6)

Recommended separation distance between portable and mobile RF communications equipment and the MDSTH1002 IR Thermometer			
MDSTH1002 IR Thermometer is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customers or the users of these MDSTH1002 IR Thermometers can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the thermometer as recommended below, according to the maximum output power of the communications equipment.			
Output Power of Transmitter in Watts	Separation distance according to frequency of transmitter, in meters		
	150 kHz to 80 MHz d = 1.2 \sqrt{P}	80 MHz to 800 MHz d = 1.2 \sqrt{P}	800 MHz to 2.5 GHz d = 2.3 \sqrt{P}
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. Note: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection for structures, objects and people.			

		<p>a. should be less than the compliance level in each frequency range.</p> <p>b. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>Note 1: At 80 MHz and 800MHz, the higher frequency range applies.</p> <p>Note 2: These guidelines may not apply to all situations.</p> <p>Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>		
<p>*Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radio, AM and FM radio broadcasts, and TV broadcasts cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the MDSTH1002 IR Thermometer is used exceeds the applicable RF compliance level above, the MDSTH1002 IR Thermometer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Clinical IR Thermometer.</p> <p>^bOver the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>		

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www.medline.com

Manufactured for: Medline Industries, Inc.,
Northfield, IL 60093 USA. Made in China
To Order Call: 1-800-MEDLINE V1 RF20DDX

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Fabricado para: Medline Industries, Inc.,