1. Remove handle (5) and slip one color-coded water jet tip (2) onto handle.
2. Remove water reservoir (3) by holding unit base firmly in place with one hand while pulling reservoir straight up with the other.
3. Fill reservoir with water, dental rinse, or mouthwash and reinstall firmly onto the unit.

Instructions for Use:

1. Prior to filling the water reservoir, make sure that the water jet tip is firmly placed onto the handle and into the stand. Once you turn on the water jet unit, it needs several seconds before the water will come out through the tip.
2. Allow full reservoir of lukewarm water to pass through hose, handle and tip before first use. (See Instructions for Use)
3. Prior to filling the water reservoir, make sure that the water jet tip is firmly placed onto the handle and into the stand. Once you turn on the water jet unit, it needs several seconds before the water will come out through the tip.

When replacing the batteries, make sure you dispose of the old batteries in accordance with applicable regulations. They should be disposed of in a location specifically designated for that purpose so they can be recycled safely and without posing any risk for the environment. Do not burn or bury them.

NOTE: Use only alkaline batteries. Do not mix old and new batteries. Do not mix alkaline, standard (carbon-zinc), or rechargeable (nickel cadmium) batteries. Remove batteries from equipment that is not to be used for an extended period of time. When replacing batteries, make sure you dispose of the old batteries properly. Clean the battery contacts and also those of the device prior to battery installation. Remove used batteries promptly.


Initial Preparation: Before using your Dental Water Jet for the first time: 1. Remove battery door (4) and insert 3 AA alkaline batteries. Ensure that the batteries are installed correctly with regard to polarity (+ and -). Replace door. Always purchase the correct size and grade of battery most suitable for intended use. Replace all batteries of a set at the same time.

2. Remove handle (5) and slip one color-coded water jet tip (2) onto handle.

NOTE: Thoroughly rinse reservoir with clean water after each use with mouthwash or dental rinse, and whenever unit has not been used for an extended period of time.

4. Begin with control at lowest setting to allow gums to become accustomed to water jet cleaning.
5. Press the PAUSE button (6), switch unit to ON, place water jet tip in mouth, and lean over the sink.
6. Release PAUSE button to begin flow and adjust HI/LO control (7) until comfortable. It usually takes several seconds before the water comes out of the water jet tip.
7. Hold mouth slightly open and elbow at mouth level to let water run back into the sink.
8. Turn handle to direct the water jet tip towards tooth surfaces and between teeth to rinse away trapped plaque and particles. The pulsing action of the water massages gums as it cleans the entire mouth.

NOTE: To prevent damage to delicate gum areas, do not point water jet tip directly at gums.


Remember: Oral irrigation is a supplement to, not a substitute for, your oral hygiene routine. Ask your dental health care provider to advise you on how the water jet can most effectively meet your specific oral hygiene needs.

The Interplak® by Conair® Dental Water Jet enhances oral hygiene and significantly contributes to healthy gums when used as part of a responsible oral health management program. It is particularly effective at cleaning around orthodontic appliances, fixed and removable bridge work, implants and splints, and it provides valuable gum massage with every use.

Welcome to a New Level of At-Home Dental Care

Cleaning and care:

If, in the course of use, the valve (9) on the water tank begins to leak, cleaning the valve seat is necessary. To do this, simply use the water jet tip to push out the valve body, clean the valve seat with warm soapy water, rinse well and replace the valve.

- | | |
|-------------------------------------|-------------------------------|
| 1 Unit | 6 Water PAUSE button |
| 2 Color-coded slip-on water jet tip | 7 HI/LO flow control |
| 3 Reservoir | 8 Power switch |
| 4 Battery door (on bottom) | 9 Valve seat (inside of tank) |
| 5 Handle | |

SAVE THESE INSTRUCTIONS

- Contact your dentist prior to use if you have had oral or periodontal surgery in the previous twelve months.
- When you start any new oral hygiene program, your gums may initially bleed slightly. Contact your dentist if excessive bleeding occurs or if bleeding continues after two weeks of use.
- Be sure that only the water jet tip touches your teeth and gums. Damage to teeth and gums may occur if any other part of the unit touches them.
- This product is designed for oral cleaning only. DO NOT use the unit for any other purpose.
- As with any other new oral hygiene program, individuals with heart problems should check with their physician prior to using the Interplak® Dental Water Jet.
- The Interplak® by Conair® Dental Water Jet is a personal care unit and is intended for household use only.

MEDICAL WARNING
For your safety:

Limited 1-Year Warranty

Conair will repair or replace (at our option) your unit free of charge for 12 months from the date of purchase if the appliance is defective in workmanship or materials.

To obtain service under this warranty, return the defective product to the service center listed below, together with your purchase receipt and \$3.00 for postage and handling.

California residents need only provide proof of purchase and should call 1-800-334-4031 for shipping instructions. In the absence of a purchase receipt, the warranty period shall be 12 months from the date of manufacture.

ANY IMPLIED WARRANTIES, OBLIGATIONS, OR LIABILITIES, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, SHALL BE LIMITED IN DURATION TO THE 12-MONTH DURATION OF THIS WRITTEN, LIMITED WARRANTY. Some states do not allow limitations on how long an implied warranty lasts, so the above limitation may not apply to you.

IN NO EVENT SHALL CONAIR BE LIABLE FOR ANY SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES FOR BREACH OF THIS OR ANY OTHER WARRANTY, EXPRESS OR IMPLIED, WHATSOEVER. Some states do not allow the exclusion or limitation of special, incidental, or consequential damages, so the above limitation may not apply to you.

This warranty gives you specific legal rights, and you may also have other rights, which vary from state to state.

Please register this product at: www.conair.com/registration

Service Center:
Conair Corporation, Service Department, 7475 N. Glen Harbor Blvd., Glendale, AZ 85307

©2012 Conair Corporation, 150 Milford Road, East Windsor, NJ 08520
For more information regarding this and other Conair products, visit www.interplak.com or call 1-800-334-4031

12IN128399
IB-11128A

SAVE THESE INSTRUCTIONS

- **DO NOT** submerge the unit in water or any other liquid.
- Close supervision is necessary when this product is used by, on, or near children or disabled persons.
- Use this product only for its intended use as described in this manual. Do not use attachments unless recommended by the manufacturer.

WARNING

DANGER

IMPORTANT SAFEGUARDS READ ALL INSTRUCTIONS BEFORE USING

INTERPLAK®

by CONAIR®

COMPACT PORTABLE WATER JET

Instruction Manual • Models WJ3CS/WJ3CSR

For your safety and continued enjoyment of this product, always read the instruction book carefully before using.

Questions or Comments? 1-800-334-4031
Visit us at our website: www.interplak.com

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Conair Corporation, Service Department, 7475 N. Glen Harbor Blvd. Glendale, AZ 85307