



Closed Wound Drainage Systems



Medline manufactures and supplies a full line of closed wound drainage products for various surgical procedures including orthopedic, plastics, general surgery, neurological and mastectomy and breast augmentation procedures.

Key products include EvaX Style Channel Silicone Wound Drains, Traditional Style Perforated Silicone Wound Drains, Silicone Bulb Evacuators, PVC Evacuator Kits and Wound Drain Kits and Penrose Drains.

EvaX Style (Channel) Silicone Wound Drains

- Utilize internal channels, or flutes, to swiftly remove exudates
- Independent channels facilitate drainage and lower risk of occlusion
- Maximum patient comfort upon removal
- Indicator dot makes drain placement quick and easy
- Available in various sizes of Round (Full-Fluted) and Flat (3/4 or Full-Fluted)
- Radiopaque

Item #	Description	Pkg
--------	-------------	-----

Flat Drains, Channel

DYNJWE2210	7mm Flat, 3/4 Fluted	10/cs
DYNJWE2211	7mm Flat, Full-Fluted	10/cs
DYNJWE2213	10mm Flat, 3/4 Fluted	10/cs
DYNJWE2214	10mm Flat, Full-Fluted	10/cs

Flat Drains w/Trocar, Channel

DYNJWE2212	7mm Flat, Full-Fluted	10/cs
DYNJWE2215	10mm Flat, Full-Fluted	10/cs
DYNJWE2216	7mm Flat, 3/4 Fluted	10/cs
DYNJWE2217	10mm Flat, 3/4 Fluted	10/cs

Flat Drain Kits, Channel

DYNJWE2205	100cc Evac w/ 7mm Full-Fluted	10/cs
DYNJWE2207	100cc Evac w/ 10mm Full-Fluted	10/cs
DYNJWE2236	100cc Evac w/7mm Flat 3/4 Fluted	10/cs
DYNJWE2238	100cc Evac w/10mm Flat 3/4 Fluted	10/cs

Round Drains, Channel


DYNJWE2186	10Fr, 1/8", Full-Fluted	10/cs
DYNJWE2188	15Fr, 3/16", Full-Fluted	10/cs
DYNJWE2190	19Fr, 1/4", Full-Fluted	10/cs
DYNJWE2234	24Fr, 5/16", Full-Fluted	10/cs

Round Drains w/Trocar, Channel

DYNJWE2187	10Fr, 1/8", Full-Fluted	10/cs
DYNJWE2189	15Fr, 3/16", Full-Fluted	10/cs
DYNJWE2191	19Fr, 1/4", Full-Fluted	10/cs

Round Drain Kits w/Trocar, Channel

DYNJWE2263	100cc Evac, 10Fr, 1/8", Full-Fluted	10/cs
DYNJWE2264	100cc Evac, 15Fr, 3/16", Full-Fluted	10/cs
DYNJWE2265	100cc Evac, 19Fr, 1/4", Full-Fluted	10/cs



Traditional Style (Perforated) Silicone Wound Drains

- Smooth perforations are heat polished for easy, painless removal.
- Perforations are positioned in alternate 90° planes for effective exudate removal
- Available in various sizes of Round and Flat drains
- Indicator dot makes drain placement quick and easy
- Radiopaque

Item #	Description	Pkg
--------	-------------	-----

Flat Drains, Perforated

DYNJWE1308	7mm x 20cm 3/4 Perforation	10/cs
DYNJWE1309	10mm x 20cm 3/4 Perforation	10/cs
DYNJWE1310	7mm x 20cm Full-perforation	10/cs
DYNJWE1311	10mm x 20cm Full-perforation	10/cs

Flat Drains w/Trocar, Perforated

DYNJWE1408	7mm x 20cm 3/4 Perforation	10/cs
DYNJWE1409	10mm x 20cm 3/4 Perforation	10/cs
DYNJWE1410	7mm x 20cm Full-perforation	10/cs
DYNJWE1411	10mm x 20cm Full-perforation	10/cs

Flat Drain Kits, Perforated

DYNJWE1348	100cc Evac w/7mm 3/4 Perforation	10/cs
DYNJWE1349	100cc Evac w/10mm 3/4 Perforation	10/cs
DYNJWE1360	100cc Evac w/7mm Full-perforation	10/cs
DYNJWE1361	100cc Evac w/10mm Full-perforation	10/cs

Flat Drain Kits w/Trocar, Perforated

DYNJWE1468	100cc Evac w/7mm 3/4 Perforation	10/cs
DYNJWE1469	100cc Evac w/10mm 3/4 Perforation	10/cs
DYNJWE1480	100cc Evac w/7mm Full Perforation	10/cs
DYNJWE1481	100cc Evac w/10mm Full Perforation	10/cs

Round Drains, Perforated

DYNJWE1320	7Fr, 3/32" End-perforation	10/cs
DYNJWE1321	10Fr, 1/8" End-perforation	10/cs
DYNJWE1323	15Fr, 3/16" End-perforation	10/cs
DYNJWE1325	19Fr, 1/4" End-perforation	10/cs

Round Drains w/Trocar, Perforated

DYNJWE0321	10Fr, 1/8" End-perforation	10/cs
DYNJWE0323	15Fr, 3/16" End-perforation	10/cs
DYNJWE0325	19Fr, 1/4" End-perforation	10/cs

Round Drain Kits w/Trocar

DYNJWE0421	100cc Evac w/10Fr, 1/8" Round	10/cs
DYNJWE0423	100cc Evac w/15Fr, 3/16" Round	10/cs
DYNJWE0425	100cc Evac w/19Fr, 1/4" Round	10/cs
DYNJWE1080	100cc Evac w/10Fr, 1/8" Round	10/cs
DYNJWE1060	100cc Evac w/15Fr, 3/16" Round	10/cs
DYNJWE1010	100cc Evac w/19Fr, 1/4" Round	10/cs



Silicone Bulb Evacuators

Designed specifically for smaller procedures where less exudate is anticipated.

Item #	Description	Pkg
DYNJWE1305	100cc capacity	10/cs
DYNJWE1000	400cc capacity	10/cs



DYNJWE1000

PVC Evacuator Kits and Wound Drains

Stronger than silicone to maximize collapse-resistance

- Available in kits, which include 400cc 3-spring evacuator, silicone tubing and "Y" connector, and drain with trocar
- Drains also available individually
- Radiopaque
- Latex-Free


Item #	Description	Pkg
400cc 3-Spring Evacuator w/PVC Drain and Trocar		
DYNJWE401	400cc w/7Fr, 3/32" PVC	12/cs
DYNJWE402	400cc w/10Fr, 1/8" PVC	12/cs
DYNJWE403	400cc w/15Fr, 3/16" PVC	12/cs
DYNJWE404	400cc w/19Fr, 1/4" PVC	12/cs

400cc 3-Spring Evacuator w/Silicone Tube Adapters

DYNJWE475	400cc 3-Spring Evac w/Y Connector	12/cs
-----------	-----------------------------------	-------

Round PVC Drains, Mid-Perforation w/Trocar

DYNJWE0521	10Fr, 1/8" Mid-perforation	10/cs
DYNJWE0523	15Fr, 3/16" Mid-perforation	10/cs
DYNJWE0525	19Fr, 1/4" Mid-perforation	10/cs



DYNJWE475

Latex Penrose Drains

Latex Penrose Drains

Item #	Description	Pkg
DYND50420	18" x 1/4" Sterile	25/cs
DYND50421	18" x 3/8" Sterile	25/cs
DYND50422	18" x 1/2" Sterile	25/cs
DYND50423	18" x 5/8" Sterile	25/cs
DYND50424	18" x 3/4" Sterile	25/cs
DYND50425	18" x 7/8" Sterile	25/cs
DYND50426	18" x 1" Sterile	25/cs

Wall Suction Adapter

For Use With 3-Spring Evacuators

Item #	Description	Pkg
DYNJWE4605	Wall Suction Adaptor	25/cs



Competitive Cross-Reference

Description	Medline	Cardinal	J&J	Zimmer	Bard
EVAX STYLE (CHANNEL) SILICONE WOUND DRAINS					
Flat Drains, Channel					
7mm Flat, 3/4 Fluted	DYNJWE2210	JP-HUF073	2210	N/A	0072210
7mm Flat, Full-Fluted	DYNJWE2211	JP-HUF070	2211	N/A	0072211
10mm Flat, 3/4 Fluted	DYNJWE2213	JP-HUF103	2213	N/A	0072213
10mm Flat, Full-Fluted	DYNJWE2214	JP-HUF100	2214	N/A	0072214
Flat Drains w/Trocar, Channel					
7mm Flat, Full-Fluted w/Trocar	DYNJWE2212	JP-HUF071	2212	N/A	0072212
10mm Flat, Full-Fluted w/Trocar	DYNJWE2215	JP-HUF101	2215	N/A	0072215
7mm Flat, 3/4 Fluted w/Trocar	DYNJWE2216	JP-HUF074	2216	N/A	0072216
10mm Flat, 3/4 Fluted w/Trocar	DYNJWE2217	JP-HUF104	2217	N/A	0072217
Flat Drain Kits, Channel					
100cc Evac w/ 7mm Full-Fluted	DYNJWE2205	JP1HGF070	2205	N/A	N/A
100cc Evac w/ 10mm Full-Fluted	DYNJWE2207	JP1HGF100	2207	N/A	N/A
100cc Evac w/7mm Flat 3/4 Fluted	DYNJWE2236	JP1HGF073	2236	N/A	N/A
100cc Evac w/10mm Flat 3/4 Fluted	DYNJWE2238	JP1HGF103	2238	N/A	N/A
Round Drains, Channel					
10fr Round, Full-Fluted	DYNJWE2186	JP-HUR100	2186 / 2226	00-2540-950-10	0072186
15fr Round, Full-Fluted	DYNJWE2188	JP-HUR150	2188 / 2228	00-2540-951-10	0072188
19fr Round Full-Fluted, 3/16" Trocar	DYNJWE2190	JP-HUR190	2190 / 2230	00-2540-952-10	0072190
19fr Round Full-Fluted, 1/4" Trocar	N/A	N/A	N/A	N/A	0072192
24fr Round Full-Fluted	DYNJWE2234	N/A	2234	N/A	N/A
Round Drains w/Trocar, Channel					
10fr Round, Full-Fluted w/Trocar	DYNJWE2187	JP-HUR101	2187 / 2221 / 2227	N/A	0072187 / 0072221
15fr Round, Full-Fluted w/Trocar	DYNJWE2189	JP-HUR151	2189 / 2223 2229 / 2233	N/A	0072189 / 0072223
19fr Round, Full-Fluted w/Trocar	DYNJWE2191	JP-HUR195 2231 / 2232	2191 / 2225	N/A	0072191 / 0072225
24fr Round, Full-Fluted w/Trocar	N/A	N/A	2234		
Round Drain Kits w/Trocar, Channel					
100cc Evacuator, 10fr Round Full-Fluted w/Trocar	DYNJWE2263	JP1HGR101	2259 / 2263	N/A	N/A
100cc Evacuator, 15fr Round Full-Fluted w/Trocar	DYNJWE2264	JP1HGR151	2264 / 2261 / 2268	N/A	N/A
100cc Evacuator, 19fr Round Full-Fluted w/Trocar	DYNJWE2265	JP1HGR195	2265 / 2262	N/A	N/A



Competitive Cross-Reference

Description	Medline	Cardinal	J&J	Zimmer	Bard
-------------	---------	----------	-----	--------	------

TRADITIONAL STYLE (PERFORATED) SILICONE WOUND DRAINS

Flat Drains, Perforated

Flat Drain 7mm x 20cm 3/4 Perforated	DYNJWE1308A	SU130-1308	N/A	00-2571-000-10	0070360
Flat Drain 10mm x 20cm 3/4 Perforated	DYNJWE1309	SU130-1309	N/A	00-2572-000-10	0070370
Flat Drain 7mm x 20cm Full Perforated	DYNJWE1310	SU130-1310	N/A	00-2569-000-10	0070430
Flat Drain 10mm x 20cm Full Perforated	DYNJWE1311	SU130-1311	N/A	00-2570-000-10	0070440

Flat Drains w/Trocar, Perforated

Flat Drain 7mm x 20cm 3/4 Perforated w/Trocar	DYNJWE1408	SU130-1408	N/A	00-2540-071-10	N/A
Flat Drain 10mm x 20cm 3/4 Perforated w/Trocar	DYNJWE1409	SU130-1409	N/A	00-2540-072-10	N/A
Flat Drain 7mm x 20cm Full w/Trocar	DYNJWE1410	SU130-1410	N/A	00-2540-069-10	N/A
Flat Drain 10mm x 20cm Full w/Trocar	DYNJWE1411	SU130-1411	N/A	00-2540-070-10	N/A

Flat Drain Kits, Perforated

100cc Evacuator w/7mm Flat 3/4 Perforated	DYNJWE1348	SU130-1348	N/A	00-2566-000-10	0071360
100cc Evacuator w/10mm Flat 3/4 Perforated	DYNJWE1349	SU130-1349	N/A	00-2567-000-10	0071370
100cc Evacuator w/7mm Flat Full Perforated	DYNJWE1360	SU130-1360	N/A	00-2564-000-10	0071430
100cc Evacuator w/10mm Flat Full Perforated	DYNJWE1361	SU130-1361	N/A	00-2565-000-10	0071440

Flat Drain Kits w/Trocar, Perforated

100cc Evacuator w/7mm 3/4 Flat w/Trocar	DYNJWE1468	N/A	N/A	00-2562-000-10	0072360
100cc Evacuator w/10mm 3/4 Flat w/Trocar	DYNJWE1469	N/A	N/A	00-2563-000-10	0072370
100cc Evacuator w/7mm Full Flat w/Trocar	DYNJWE1480	N/A	N/A	00-2560-000-10	0072430
100cc Evacuator w/10mm Full Flat w/Trocar	DYNJWE1481	N/A	N/A	00-2561-000-10	0072440

Round Drains, Perforated

Round Drain 3/32" Endperforated	DYNJWE1320	SU130-1320	N/A	00-2540-001-10	0070200
Round Drain 1/8" Endperforated	DYNJWE1321A	SU130-1321	N/A	00-2540-002-10	0070210
Round Drain 3/16" Endperforated	DYNJWE1323A	SU130-1323	N/A	00-2540-007-10	0070220
Round Drain 1/4" Endperforated	DYNJWE1325A	SU130-1325	N/A	00-2540-008-10	0070230

Round Drains w/Trocar, Perforated

Round Drain 1/8" Endperforated w/Trocar 10" Hole Pattern	DYNJWE0321	SU130-0321	N/A	00-2540-022-10	0070310
Round Drain 3/16" Endperforated w/Trocar, 10" Hole Pattern	DYNJWE0323	SU130-0323	N/A	00-2540-023-10	0070320
Round Drain 1/4" Endperforated w/Trocar 10" Hole Pattern	DYNJWE0325	SU130-0325	N/A	00-2540-024-10	0070330

Round Drain Kits with Trocar, Perforated

100cc Evacuator w/ 1/8" Round and Trocar	DYNJWE0421	N/A	N/A	00-2536-000-10	0073310
100cc Evacuator w/ 3/16" Round and Trocar	DYNJWE0423	N/A	N/A	00-2537-000-10	0073320
100cc Evacuator w/ 1/4" Round and Trocar	DYNJWE0425	N/A	N/A	00-2538-000-10	0073330
400cc Evacuator w/ 1/8" Round and Trocar	DYNJWE1080	SU130-1080	N/A	N/A	N/A
400cc Evacuator w/ 3/16" Round and Trocar	DYNJWE1060	SU130-1060	N/A	N/A	N/A
400cc Evacuator w/ 1/4" Round and Trocar	DYNJWE1010	SU130-1010	N/A	N/A	N/A

Competitive Cross-Reference

Description	Medline	Cardinal	J&J	Zimmer	Bard
SILICONE BULB EVACUATORS					
Silicone Bulb Evacuator					
100cc Silicone Evacuator	DYNJWE1305	SU130-1305	2160	00-2568-00-10	0070740
400cc Silicone Evacuator	DYNJWE1000	SU130-1000	N/A	N/A	N/A

PVC EVACUATOR KITS AND WOUND DRAINS					
400cc 3-Spring Evacuator w/Silicone Tube Adapters					
400cc 3-Spring Evacuator w/Y Connector	DYNJWE475	SU130-475	N/A	00-2500-004-10	0043650
400cc 3-Spring Evacuator w/Y Connector	DYNJWE475	SU130-475	N/A	00-2500-004-10	0070800

Round PVC Drains, Mid-Perforation w/Trocar					
Pvc Round Drain 1/8" Midperforation w/5" Trocar	DYNJWE0521	SU130-0521	N/A	00-1500-022-10	0080710
Pvc Round Firm Drain 3/16" Midperforation w/5" Trocar	DYNJWE0523	SU130-0523	N/A	00-1500-024-10	0080720
Pvc Round Drain 1/4" Midperforation w/5" Trocar	DYNJWE0525	SU130-0525	N/A	00-1500-025-10	0080730

400cc 3-Spring Evacuator w/PVC Drain and Trocar					
400cc 3-Spring Evacuator w/ 3/32" PVC Trocar	DYNJWE401	SU130-401D	N/A	00-2595-000-10	0043600
400cc 3-Spring Evacuator w/ 1/8" PVC Drain and Trocar	DYNJWE402	SU130-402D	N/A	00-2500-000-10	0043610
400cc 3-Spring Evacuator w/ 3/16" PVC Drain and Trocar	DYNJWE403	SU130-403D	N/A	00-2505-000-10	0043630
400cc 3-Spring Evacuator w/ 1/4" PVC Drain and Trocar	DYNJWE404	SU130-404D	N/A	00-2590-000-10	0043620
400cc 3-Spring Evacuator w/ 1/4" PVC Drain and 5" Trocar	DYNJWE404	N/A	N/A	00-2590-040-10	N/A

WALL SUCTION ADAPTOR: FOR USE W/3-SPRING EVACUATORS					
Wall Suction Adaptor	DYNJWE4605	SU130-4605	N/A	00-1500-060-10	0034850

PENROSE DRAINS					
Latex Penrose Drain 18" x 1/4" Sterile	DYND50420	20416-025	N/A	N/A	0918010
Latex Penrose Drain 18" x 3/8" Sterile	DYND50421	N/A	N/A	N/A	0918020
Latex Penrose Drain 18" x 1/2" Sterile	DYND50422	20416-050	N/A	N/A	0918030
Latex Penrose Drain 18" x 5/8" Sterile	DYND50423	N/A	N/A	N/A	0918040
Latex Penrose Drain 18" x 3/4" Sterile	DYND50424	20416-075	N/A	N/A	0918050
Latex Penrose Drain 18" x 7/8" Sterile	DYND50425	N/A	N/A	N/A	0918060
Latex Penrose Drain 18" x 1" Sterile	DYND50426	20416-100	N/A	N/A	0918070



www.medline.com 1-800-MEDLINE (1-800-633-5463)

©2004 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc. One Medline Place, Mundelein, IL 60060
MKT204050/Lit127/5M