

AAMI Standards for Proxima Surgical Gowns

LEVEL 1 *No gowns available*

LEVEL 2 **Ability to resist liquid penetration in impact penetration and hydrostatic pressure tests in critical zone areas.**

Eclipse® Nonreinforced

DYNJP2001	Large
DYNJP2002	X-Large
DYNJP2003	XX-Large

LEVEL 3 **Ability to resist liquid penetration in impact penetration and hydrostatic pressure tests in critical zone areas. (hydrostatic pressure at higher level)**

Eclipse® Fabric Reinforced

DYNJP2101	Large
DYNJP2102	X-Large
DYNJP2103	XX-Large

Sirus™ Nonreinforced

DYNJP2401	Large
DYNJP2402	X-Large
DYNJP2403	XX-Large

Sirus™ Fabric Reinforced

DYNJP2501	Large
DYNJP2502	X-Large
DYNJP2503	XX-Large

Aurora™ Nonreinforced

DYNJP2701	Large
DYNJP2702	X-Large
DYNJP2703	XX-Large

Aurora™ Nonreinforced, Raglan Sleeve

DYNJP2801	Large
DYNJP2802	X-Large
DYNJP2803	XX-Large

Aurora™ Fabric Reinforced

DYNJP2704	Large
DYNJP2705	X-Large
DYNJP2706	XX-Large

Aurora™ Fabric Reinforced, Raglan Sleeve

DYNJP2804	Large
DYNJP2805	X-Large
DYNJP2806	XX-Large

LEVEL 4 **Ability to resist liquid and viral penetration in all critical zone areas (ASTM F1671)**

Eclipse® Poly Reinforced

DYNJP2201	Large
DYNJP2202	X-Large
DYNJP2203	XX-Large

Eclipse® Poly Reinforced, X-Long

DYNJP2204	Large
DYNJP2205	X-Large
DYNJP2206	XX-Large

Eclipse® Poly Reinforced with Breathable, Impervious Sleeves

DYNJP2207	Large
DYNJP2208	X-Large
DYNJP2209	XX-Large

Sirus™ Poly Reinforced

DYNJP2601	Large
DYNJP2602	X-Large
DYNJP2603	XX-Large

Aurora™ Poly Reinforced

DYNJP2707	Large
DYNJP2708	X-Large
DYNJP2709	XX-Large

Aurora™ Poly Reinforced, Raglan Sleeve

DYNJP2807	Large
DYNJP2808	X-Large
DYNJP2809	XX-Large

Prevention™ Breathable Impervious

DYNJP2301	Large
DYNJP2302	X-Large, X-Long
DYNJP2303	XX-Large, X-Long

