



TRADITIONAL WOUND CARE

Reference Guide



VALUE AND QUALITY
FOR HEALTH CARE



COTTON WOVEN GAUZE SPONGES



STERILE

Ideal for wound debriding, prepping, packing, dressing and general wound care.

- 100% cotton sponges constructed of fine mesh gauze
- Low lint manufacturing process assures debris-free packing
- "C fold" design minimizes loose threads and lint

Sterile 2's

Item Number	Item Desc.	Case Pkg
PRM2208	2" x 2", 8-ply	50 pk/bx, 30 bx/cs, 1,500 pks/cs
PRM2208Z	2" x 2", 8-ply	50 pk/bx
PRM3312	3" x 3", 12-ply	40 pk/bx, 30 bx/cs, 1,200 pks/cs
PRM3312Z	3" x 3", 12-ply	40 pk/bx
PRM4408	4" x 4", 8-ply	50 pk/bx, 12 bx/cs, 600 pks/cs
PRM4408Z	4" x 4", 8-ply	50 pk/bx
PRM9408/Econ	4" x 4", 8 ply	100/bx, 12 bx/cs, 1,200/cs
PRM4412	4" x 4", 12-ply	25 pk/bx, 24 bx/cs, 600 pks/cs
PRM4412H	4" x 4", 12-ply	25 pk/bx
PRM9412/Econ	4" x 4", 12-ply	25 pk/bx, 24 bx/cs, 600 pks/cs

Sterile 1's

Item Number	Item Desc.	Case Pkg
PRM21419	2" x 2", 12-ply	100 pk/bx, 24 bx/cs, 2,400 pks/cs
PRM21423	4" x 4", 12-ply	100 pk/bx, 12 bx/cs, 1,200 pks/cs

Q: When should I use cotton gauze?

- ✓ Cleaning
- ✓ Debriding
- ✓ Packing
- ✓ Scrubbing
- ✓ Prepping

NON-STERILE

The most cost-effective choice for cleaning, debriding and covering wounds.

- Packaged in sleeves of 200 sponges
- The preferred choice of most nursing homes

Item Number	Item Description	Case Pkg
PRM21208C	2" x 2", 8-ply	200/sleeve, 25 pk/cs, 5,000/cs
PRM21208CZ	2" x 2", 8-ply	200/sleeve
PRM21212C	2" x 2", 12-ply	200/sleeve, 40 pk/cs, 8,000/cs
PRM21312C	3" x 3", 12-ply	200/sleeve, 20 pk/cs, 4,000/cs
PRM21312CZ	3" x 3", 12-ply	200/sleeve
PRM21408C	4" x 4", 8-ply	200/sleeve, 20 pk/cs, 4,000/cs
PRM21408CZ	4" x 4", 8-ply	200/sleeve
PRM91408C	4" x 4", 8-ply, economy	200/sleeve, 20 pk/cs, 4,000/cs
PRM21412C	4" x 4", 12-ply	200/sleeve, 10 pk/cs, 2,000/cs
PRM91412C	4" x 4", 12-ply economy	200/sleeve, 10 pk/cs, 2,000/cs
PRM21416C	4" x 4", 16-ply	200/sleeve, 10 pk/cs, 2,000/cs

Cotton Filled Dental Sponges

PRM22DS	2" x 2"	200/sleeve, 25 pk/cs, 5,000/cs
---------	---------	--------------------------------



NON-WOVEN GAUZE SPONGES



More absorbency and less linting than 100% cotton gauze.

- Unique rayon/polyester formed fabric
- Significantly less linting than cotton gauze
- Virtually non-adherent
- More cost-effective than cotton gauze

Sterile 2's

Item Number	Item Description	Case Pkg
PRM21224	2" x 2", 4-ply, 2's	50pk/bx, 30 bx/cs, 1500 pks/cs
PRM21444	4" x 4", 4-ply, 2's	50pk/bx, 12 bx/cs, 600 pks/cs

Non-Sterile - Bulk

Item Number	Item Desc.	Case Pkg
PRM25224	2" x 2", 4-ply	200/sleeve, 40 sleeves/cs, 8000/cs
PRM25444	4" x 4", 4-ply	200/sleeve, 10 sleeves/cs, 2000/cs

Sterile 1's

Item Number	Item Description	Case Pkg
DG21334	3" x 3", 4-ply	1ea/pouch, 10 pouches/box, 24ea/cs

Q: When should I use non-woven?

- ✓ Absorption
- ✓ Low-linting
- ✓ Fast wicking
- ✓ Wound release

ABD/COMBINE PADS

For heavily draining wounds.

- Soft, non-woven rayon/polyester outer layer is virtually non-adherent
- Compressed cellulose center quickly absorbs heavy exudate
- Hydrophobic back retards strike through
- Sealed, finely finished edges prevent linting

STERILE 1's

Item Number	Item Desc.	Case Pkg
PRM21450	5" x 9"	25/bx, 16 bx/cs, 400/cs
PRM21450Z	5" x 9"	25/bx
PRM21453	8" x 7.5"	20/bx, 12 bx/cs, 240/cs
PRM21454	8" x 10"	16/bx, 20 bx/cs, 320/cs

NON-STERILE Bulk

Item Number	Item Desc.	Case Pkg
PRM21451	5" x 9"	576/cs



Most Carling Traditional Wound Care products are latex-free

To order, call 1-800-633-5463

COTTON GAUZE BANDAGE ROLLS

Ideal for bandaging heads, limbs and difficult to dress wounds.

- 6-ply, 100% cotton
- Unique crimped, open-weave design means superior absorbency and cushioning
- Gently stretches and conforms to body contours
- Used as either a primary or secondary dressing

Item Number	Item Desc.	Case Pkg
PRM25865	Sterile, 4.5" x 4.1 yds	100 ea/cs
PRM25865H	Sterile, 4.5" x 4.1yds	1 ea
PRM25855	Non-sterile, 4.5" x 4.1 yds	100 ea/cs
DG25865	Sterile, 4.5" x 4 yds	1ea/pouch, 24ea/cs



SUPRA FORM™ CONFORMING BANDAGES



Gently secures primary dressings without restriction of movement.

- One-ply, rayon/polyester weave conforms to hard to wrap areas
- One-way stretch with memory accommodates swelling
- Finished edges for low linting

STERILE - Unstretched

Item Number	Item Desc.	Case Pkg
PRM25496	2" x 75"	12/bx, 8 bx/cs, 96/cs
PRM25497	3" x 75"	12/bx, 8 bx/cs, 96/cs
PRM25497Z	3" x 75"	12/bx
PRM25498	4" x 75"	12/bx, 8 bx/cs, 96/cs
PRM25498Z	4" x 75"	12/bx
PRM25499	6" x 80"	8/bx, 6 bx/cs, 48/cs

NON-STERILE - Unstretched

Item Number	Item Desc.	Case Pkg
PRM25492	2" x 75"	12/bag, 8 bag/cs, 96/cs
PRM25493	3" x 75"	12/bag, 8 bags/cs, 96/cs
PRM25493Z	3" x 75"	12/bag
PRM25494	4" x 75"	12/bag, 8 bags/cs, 96/cs
PRM25494Z	4" x 75"	12/bag
PRM25495	6" x 80"	12/bag, 4 bags/cs, 48/cs


DRAIN SPONGES

Unique opening helps to wick away more fluid from tube site.

- Superior fenestration more closely hugs drain tubes, keeping tube sites drier
- 6-ply, rayon/polyester fabric provides maximum absorbency and low linting
- Pre-cut sponge eliminates the need for sterile scissors and conserves nursing time

Sterile 2's

Item Number	Item Desc.	Case Pkg
PRM256001	2" x 2", 6-ply	35 pk/bx, 20 bx/cs, 700 pks/cs
PRM256000	4" x 4", 6-ply	25 pk/bx, 12 bx/cs, 300 pks/cs
PRM256000H	4" x 4", 6-ply	25 pk/bx



CARING SELF-ADHERENT WRAP



- Self-adherent, (cohesive) wrap, sticks to itself but not to skin, hair or dressings. No pins or clips required.
- Secondary dressing that functions like a tape to retain dressings without tension around digits or limbs.
- Can be applied to provide light compression if desired.
- Pain free removal that will not damage fragile skin.
- Comfortable, porous and lightweight construction encourages patient compliance.
- Ideal for wrapping splints, gauze pads and dressings securement, compression wrap post blood draw, I.V. securement
- Each bandage is individually wrapped.

Tan, LATEX-FREE

Item Number	Item Desc.	Case Pkg
PRM088001	1" X 5 yds stretched, Latex-free	30 rolls/cs
PRM088002	2" X 5 yds stretched, Latex-free	36 rolls/cs
PRM088003	3" X 5 yds stretched, Latex-free	24 rolls/cs
PRM088004	4" X 5 yds stretched, Latex-free	18 rolls/cs
PRM088006	6" X 5 yds stretched, Latex-free	12 rolls/cs

Tan, LATEX

Item Number	Item Desc.	Case Pkg
PRM086001	1" X 5 yds stretched, Latex	30 rolls/cs
PRM086002	2" X 5 yds stretched, Latex	36 rolls/cs
PRM086003	3" X 5 yds stretched, Latex	24 rolls/cs
PRM086004	4" X 5 yds stretched, Latex	18 rolls/cs
PRM086006	6" X 5 yds stretched, Latex	12 rolls/cs

LARGE ADHESIVE PAD



- Sheer plastic backing is lightweight, conformable, breathable and designed to be less noticeable on skin.
- Non-adherent, absorbent, island pad protects granulating tissue.
- Durable adhesion seals out dirt.
- Apply to clean, dry skin.
- Individually wrapped.
- Latex free
- Sterile

Item Number	Item Desc.	Case Pkg
DG0726	3" X 4"	10 ea/box, 24bx/case, 240ea/cs

To order, call 1-800-633-5463

PLASTIC ADHESIVE BANDAGES



Protect minor wounds and scrapes.

- Tan plastic backing is water resistant
- Long-lasting adhesive with minimal residue
- Absorbent non-adherent continuous pad
- Sterile

Item Number	Item Desc.	Case Pkg
PRM25500	3/4" x 3"	100/bx, 12 bx/cs, 1200/cs
PRM25500H	3/4" x 3"	100/bx
PRM25600	1" x 3"	100/bx, 12 bx/cs, 1200/cs
PRM25600H	1" x 3"	100/bx

CARING SURGICAL TAPES

TRANSPARENT TAPE

Permits skin examination without tape removal.

- Clear, breathable plastic material
- Micro-perforation allows for bi-directional tearing
- Excellent for IV sites and tubing retention

Item Number	Item Desc.	Case Pkg
PRM260201	1" x 10 yds	12 rolls/bx, 12 bx/cs, 144 rolls/cs
PRM260201H	1" x 10 yds	12 rolls/bx
PRM260202	2" x 10 yds	6 rolls/bx, 12 bx/cs, 72 rolls/cs
PRM260202H	2" x 10 yds	6 rolls/bx
PRM260203	3" x 10 yds	4 rolls/bx, 12 bx/cs, 48 rolls/cs
PRM260203H	3" x 10 yds	4 rolls/bx

CLOTH "SILK" TAPE

A strong tape for affixing dressings.

- Sturdy rayon/acetate fabric helps secure dressings in place
- Serrated edges facilitate tearing
- Durable adhesion with minimal residue on skin
- Water repellent

Item Number	Item Desc.	Case Pkg
PRM260101	1" x 10 yds	12 rolls/bx, 12 bx/cs, 144 rolls/cs
PRM260101H	1" x 10 yds	12 rolls/bx
PRM260102	2" x 10 yds	6 rolls/bx, 12 bx/cs, 72 rolls/cs
PRM260102H	2" x 10 yds	6 rolls/bx
PRM260103	3" x 10 yds	4 rolls/bx, 12 bx/cs, 48 rolls/cs



PAPER TAPE

Gentle, breathable tape for sensitive skin.

- Gentle, air-permeable, non-woven material
- Especially recommended for geriatric patients or others with fragile skin
- Conforms and adheres even to damp skin
- Secures lightweight dressing and tubing, ostomy pouches and facial dressings

Item Number	Item Desc.	Case Pkg
PRM260001	1" x 10 yds	12 rolls/bx, 12 bx/cs, 144 rolls/cs
PRM260001H	1" x 10 yds	12 rolls/bx
PRM260002	2" x 10 yds	6 rolls/bx, 12 bx/cs, 72 rolls/cs
PRM260002H	2" x 10 yds	6 rolls/bx
PRM260003	3" x 10 yds	4 rolls/bx, 12 bx/cs, 48 rolls/cs
PRM260003H	3" x 10 yds	4 rolls/bx

Need Traditional Wound Care items you can't find in this brochure?



If you can't find something you need, please see Medline's Traditional Wound Care brochure or ask your Medline sales representative for assistance. Call 1-800-MEDLINE (633-5463), or visit our on-line catalog at www.medline.com.

1-800-MEDLINE (633-5463) www.medline.com

©2008 Medline Industries, Inc., Medline is a registered trademark of Medline Industries, Inc., One Medline place, Mundelein, IL 60060. We reserve the right to correct any errors that may occur within this brochure.

Avant®, Carling® and Sheer-Gard® are registered trademarks of Medline Industries, Inc. Supra Form™ is a trademark of Medline Industries, Inc.



Most Carling Traditional Wound Care products are latex-free

MKT207376/LIT163/10M/4ZPN