




A Commitment to Quality.

QC You Can See,
PT/INR Results You Can Trust

The INRatio2 monitor offers two levels of on board quantitative quality controls for testing confidence and accuracy, making it an excellent and convenient choice for in-office or in-home use. The INRatio2 performs normal and therapeutic assays automatically on every test and provides quantitative, numeric, recordable results on the display – true QC you can see.



See how the competition stacks up.



The INRatio2® delivers the information that you rely on to manage anticoagulation patients in your office. No call backs, no waiting—test patients and direct any dosage changes face-to-face at the point of care.

Product Features

	INRatio2 ¹	CoaguChek®XS	ITC Pro Time 3®
60-second test time	yes	yes	no ²
Recordable QC results	yes	no	no
Individually wrapped strips	yes	no	yes
200-test battery life	yes	no ³	no
120-results memory	yes	no	no
Hematocrit range	25-55%	25-55%	20-60%

You're in Good Company

The #1 Selling Point of Care PT/INR Monitor for Three Years in a Row.⁴

The INRatio2 testing system combines reliable results with practical convenience, making it an excellent in-office solution for anticoagulation management.



Medline Ordering Information

P-T200431 INRatio 2 Professional Kit (includes Meter)
 P-T100139Z INRatio Test Strips 48/box
 P-T0200235A MicroSafe Capillary Tubes 100/bag
 OMMAT0714Z Safety Lancet, 21 Gauge 200/box

Medline Industries, Inc.
 One Medline Place
 Mundelein, IL 60060
 1-800-MEDLINE (1-800-633-5463)
www.medline.com/products/ptinr/

Smaller, Faster and Easier to Use⁵

- Larger LCD screen and easy-read digits
- Simplified, icon-based interface
- Battery life supports 200 tests
- Expanded test strip area for easier sample application
- Compatible with INRatio Test Strips
- Auto turn-on with strip insertion
- 33% smaller than INRatio

Sources

¹ INRatio2 Monitoring System package insert, San Jose, CA. HemoSense, Inc. 2007

² ProTime® System package insert, Edison, NJ. ITC, 2006

³ CoaguChek®XS System package insert, Indianapolis, IN. Roche Diagnostics Corporation, 2006


⁴ HPIS comparative sales data, point of care PT/INR monitors sold Q2 2004 - Q3 2007

⁵ Compared to INRatio

PN 5500327 RevC

CoaguChek® is a registered trademark of Roche Diagnostics. ProTime® is a registered trademark of ITC.

HemoSense® and INRatio® are registered trademarks of the Inverness Medical group of companies. ©2008 Inverness Medical. All rights reserved.

an  inverness medical company