

A close-up photograph of a surgeon in an operating room. The surgeon is wearing a blue bouffant cap, clear safety goggles, and a light blue surgical mask. They are wearing a blue surgical gown and green gloves. The background is a blurred operating room with blue drapes and surgical equipment.

**COMFORT & QUALITY
THAT PERFORMS.**

Surgical Gowns
The Largest Selection in the Industry.



**DESIGNED, TESTED
AND TRUSTED
BY O.R.
PROFESSIONALS**

MEDLINE SURGICAL GOWNS

Performance, variety and quality

Comfortable performance

We've designed our surgical gowns to enhance the performance of your surgical team. They're cut wide in the chest and sleeves for greater comfort, easier donning and freedom of movement. Neckbands use hook-and-loop closures for an easily adjustable, gap-free custom fit. A wide range of sizes and lengths satisfy the needs of every person in the O.R. High-tech, breathable fabric keeps the wearer comfortable and dry, even during long procedures.

Abundant choices

You want options in style, fabric, sizing and reinforcement to satisfy any preference in the O.R. Medline's four gown families each feature their own distinct fabric. You have choices in sizes, lengths and sleeve construction. We're proud to say that over 4 million Medline gowns are sold to satisfied customers every month.

Dependable quality

Our gown fabrics offer incredible choice. Low-linting and not made with natural rubber latex. Aurora, Eclipse and Sirius Gowns meet or exceed recommendations for fluid resistance from the Association for the Advancement of Medical Instrumentation (AAMI PB70). All gown material meets or exceeds the flammability guidelines found in CPSC's 16 CFR, part 1610.

SURGICAL GOWNS

AURORA™ SURGICAL GOWNS

Unbeatable combination of AAMI Level 3 and AAMI Level 4 protection and comfort

4-5

ECLIPSE® SURGICAL GOWNS

Light and breathable, for maximum comfort

6-7

SIRUS® SURGICAL GOWNS

The industry standard

8-9

PREVENTION™ PLUS GOWNS

All-around breathable film coverage

10-11



AURORATM SURGICAL GOWNS

An unbeatable combination of AAMI Level 3 and AAMI Level 4 protection and comfort

The choice for a high level of AAMI Level 3 and AAMI Level 4 protection

We understand comfort and protection are important traits in a surgical gown, especially during extended procedures. But most gowns require a compromise of one quality or the other. Aurora surgical gowns are the exception. They combine both exceptional barrier performance and superb breathability, giving surgeons the best of both worlds.

Excellent AAMI Level 3 and AAMI Level 4 protection

Exceptional barrier performance with a high hydrostatic head.

Breathable and comfortable

Advanced material quickly transfers heat away, keeping the wearer cool and comfortable.

Optional raglan cut

Roominess and a larger sweep makes movement easier.



Non-reinforced	Size	Non-Sterile Component
DYNJP2701: 30/cs	Large (43", 109 cm)	SPT-2701K
DYNJP2702: 30/cs	X-Large (47", 119 cm)	SPT-2702K
DYNJP2703: 18/cs	XX-Large (49", 124 cm)	SPT-2703K
DYNJP2715: 30/cs	Small (41", 104 cm)	SPT-2715K



AAMI LEVEL 3

Non-reinforced, Raglan Sleeve	Size	Non-Sterile Component
DYNJP2801: 32/cs	Large (45", 114 cm)	SPT-2801K
DYNJP2802: 30/cs	X-Large (48", 122 cm)	SPT-2802K
DYNJP2803: 28/cs	XX-Large (50", 127 cm)	SPT-2803K



AAMI LEVEL 3

Fabric-reinforced	Size	Non-Sterile Component
DYNJP2704: 30/cs	Large (43", 109 cm),	SPT-2704K
DYNJP2705: 30/cs	X-Large (47", 119 cm)	SPT-2705K
DYNJP2706: 18/cs	XX-Large (49", 124 cm)	SPT-2706K



AAMI LEVEL 3

Fabric-reinforced, Raglan Sleeve	Size	Non-Sterile Component
DYNJP2804: 30/cs	Large (45", 114 cm)	SPT-2804K
DYNJP2805: 28/cs	X-Large (48", 122 cm)	SPT-2805K
DYNJP2806: 28/cs	XX-Large (50", 127 cm)	SPT-2806K



AAMI LEVEL 3

Poly-reinforced	Size	Non-Sterile Component
DYNJP2707: 30/cs	Large (43", 109 cm)	SPT-2707K
DYNJP2708: 30/cs	X-Large (47", 119 cm)	SPT-2708K
DYNJP2709: 18/cs	XX-Large (49", 124 cm)	SPT-2709K



AAMI LEVEL 4

Poly-reinforced, Extra-long	Size	Non-Sterile Component
DYNJP2724: 28/cs	Large, Extra-Long (52", 132 cm)	SPT-2724K
DYNJP2725: 28/cs	X-Large, Extra-Long (56", 142 cm)	SPT-2725K
DYNJP2726: 18/cs	XX-Large, Extra-Long (56", 142 cm)	SPT-2726K



AAMI LEVEL 4

Poly-reinforced, Raglan Sleeve	Size	Non-Sterile Component
DYNJP2807: 32/cs	Large (45", 114 cm)	SPT-2807K
DYNJP2808: 30/cs	X-Large (48", 122 cm)	SPT-2808K
DYNJP2809: 28/cs	XX-Large (50", 127 cm)	SPT-2809K



AAMI LEVEL 4

Poly-reinforced, Breathable, Impervious Sleeves	Size	Non-Sterile Component
DYNJP2727: 30/cs	Large (43", 109 cm)	SPT-2727K
DYNJP2728: 30/cs	X-Large (47", 119 cm)	SPT-2728K
DYNJP2729: 18/cs	XX-Large (49", 124 cm)	SPT-2729K
DYNJP2712: 18/cs	XXX-Large (49", 124 cm)	SPT-2712K



AAMI LEVEL 4

ECLIPSE[®] SURGICAL GOWNS

Light and breathable
for maximum comfort

Everyday comfort

Appreciate the comfort of Eclipse. Using the latest advancements in SMS technology, Eclipse gowns offer outstanding softness and enhanced breathability. Stay cool and comfortable while you work.

Advanced material technology

Latest generation, softened SMS lightweight technology provides exceptional comfort.

Extended sizing

Gowns up to size 4X. Regular and extra-long lengths.

Light weight

Great protection packaged in a lightweight material with improved wicking and moisture evaporation.



Non-reinforced	Size	Non-Sterile Component
DYNJP2001: 30/cs	Large (43", 109 cm)	SPT-2001E
DYNJP2002: 30/cs	X-Large (47", 119 cm)	SPT-2002E
DYNJP2003: 18/cs	XX-Large (49", 124 cm)	SPT-2003E
DYNJP2004: 18/cs	XXX-Large (50", 127 cm)	SPT-2004E
DYNJP2005: 30/cs	Small (41", 104 cm)	SPT-2005E
DYNJP2009: 18/cs	XXXX-Large (49", 124 cm)	SPT-2009E



AAMI LEVEL 2

Fabric-reinforced	Size	Non-Sterile Component
DYNJP2101: 30/cs	Large (43", 109 cm)	SPT-2101E
DYNJP2102: 30/cs	X-Large (47", 119 cm)	SPT-2102E
DYNJP2103: 18/cs	XX-Large (49", 124 cm)	SPT-2103E



AAMI LEVEL 3

Poly-reinforced	Size	Non-Sterile Component
DYNJP2201: 30/cs	Large (43", 109 cm)	SPT-2201E
DYNJP2202: 30/cs	X-Large (47", 119 cm)	SPT-2202E
DYNJP2203: 18/cs	XX-Large (49", 124 cm)	SPT-2203E



AAMI LEVEL 4

Poly-reinforced, Extra-long	Size	Non-Sterile Component
DYNJP2204: 28/cs	Large, Extra-Long (52", 132 cm)	SPT-2204E
DYNJP2205: 28/cs	X-Large, Extra-Long (56", 142 cm)	SPT-2205E
DYNJP2206: 18/cs	XX-Large, Extra-Long (56", 142 cm)	SPT-2206E

Poly-reinforced, Breathable Impervious Sleeve	Size	Non-Sterile Component
DYNJP2207: 30/cs	Large (43", 109 cm)	SPT-2207E
DYNJP2208: 30/cs	X-Large (47", 119 cm)	SPT-2208E
DYNJP2209: 18/cs	XX-Large (49", 124 cm)	SPT-2209E
DYNJP2212: 18/cs	XXX-Large (50", 127 cm)	SPT-2212E



AAMI LEVEL 4

SIRUS[®] SURGICAL GOWNS

The industry standard

Reliability. Comfort. Period.

With some of the most advanced SMS fabric available, Sirus surgical gowns offer a high level of comfort. Dependable AAMI Level 3 and AAMI Level 4 protection comes from five fabric layers: three densely-packed meltblown layers sandwiched between two strong spunbound outer layers, for an excellent barrier against fluids and pathogens.

Extended sizing

Up to 4XL.

Optional raglan cut

Roominess and a larger sweep make movement easier.

The most choices

With over 30 variations on size, style and protection options, you'll find the right gown for the person and procedure every time.



Non-reinforced	Size	Non-Sterile Component
DYNJP2001S: 20/cs	Large (43", 109 cm)	SPT-2001SMS
DYNJP2002S: 20/cs	X-Large (47", 119 cm)	SPT-2002SMS
DYNJP2003S: 18/cs	XX-Large (49", 124 cm)	SPT-2003SMS
DYNJP2005S: 30/cs	Small (41", 104 cm)	SPT-2005SMS
DYNJP2004S: 18/cs	XXX-Large (49", 124 cm)	SPT-2004SMS
DYNJP2009S: 18/cs	XXXX-Large (49", 124 cm)	SPT-2009SMS



AAMI LEVEL 3

Non-reinforced, Extra-long	Size	Non-Sterile Component
DYNJP2002SL: 20/cs	X-Large, Extra-Long (56", 142 cm)	SPT-2002SL
DYNJP2003SL: 18/cs	XX-Large, Extra-Long (56", 142 cm)	SPT-2003SL



AAMI LEVEL 3

Non-reinforced Raglan Sleeve	Size	Non-Sterile Component
DYNJP2401: 32/cs	Large (45", 114 cm)	SPT-2401SMS
DYNJP2402: 30/cs	X-Large (48", 122 cm)	SPT-2402SMS
DYNJP2403: 28/cs	XX-Large (50", 127 cm)	SPT-2403SMS



AAMI LEVEL 3

Fabric-reinforced	Size	Non-Sterile Component
DYNJP2101S: 20/cs	Large (43", 109 cm)	SPT-2101SMS
DYNJP2102S: 20/cs	X-Large (47", 119 cm)	SPT-2102SMS
DYNJP2103S: 18/cs	XX-Large (49", 124 cm)	SPT-2103SMS



AAMI LEVEL 3

Fabric-reinforced Raglan Sleeve	Size	Non-Sterile Component
DYNJP2501: 30/cs	Large (45", 114 cm)	SPT-2501SMS
DYNJP2502: 28/cs	X-Large (48", 122 cm)	SPT-2502SMS
DYNJP2503: 28/cs	XX-Large (50", 127 cm)	SPT-2503SMS

Poly-reinforced	Size	Non-Sterile Component
DYNJP2201S: 20/cs	Large (43", 109 cm)	SPT-2201SMS
DYNJP2202S: 20/cs	X-Large (47", 119 cm)	SPT-2202SMS
DYNJP2203S: 18/cs	XX-Large (49", 124 cm)	SPT-2203SMS



AAMI LEVEL 4

Poly-reinforced, Extra-long	Size	Non-Sterile Component
DYNJP2204S: 20/cs	Large, Extra-Long (52", 132 cm)	SPT-2204SMS
DYNJP2205S: 20/cs	X-Large, Extra-Long (56", 142 cm)	SPT-2205SMS
DYNJP2206S: 18/cs	XX-Large, Extra-Long (56", 142 cm)	SPT-2206SMS
DYNJP2219SL: 18/cs	XXX-Large, Extra-Long (56", 142 cm)	SPT-2219SL

Poly-reinforced Raglan Sleeve	Size	Non-Sterile Component
DYNJP2601: 32/cs	Large (45", 114 cm)	SPT-2601SMS
DYNJP2602: 30/cs	X-Large (48", 122 cm)	SPT-2602SMS
DYNJP2603: 28/cs	XX-Large (50", 127 cm)	SPT-2603SMS

Poly-reinforced, Breathable Impervious Sleeves	Size	Non-Sterile Component
DYNJP2207S: 30/cs	Large (43", 109 cm)	SPT-2207SMS
DYNJP2208S: 30/cs	X-Large (47", 119 cm)	SPT-2208SMS
DYNJP2209S: 18/cs	XX-Large (49", 124 cm)	SPT-2209SMS



AAMI LEVEL 4

Poly-reinforced, Breathable Impervious Sleeves, Extra-long	Size	Non-Sterile Component
DYNJP2227S: 20/cs	Large, Extra-Long (52", 132 cm)	SPT-2227SMS
DYNJP2228S: 20/cs	X-Large, Extra-Long (56", 142 cm)	SPT-2228SMS
DYNJP2229S: 18/cs	XX-Large, Extra-Long (56", 142 cm)	SPT-2229SMS
DYNJP2212SL: 18/cs	XXX-Large, Extra-Long (56", 142 cm)	SPT-2212SL



AAMI LEVEL 4

PREVENTION™ PLUS GOWNS

All-around breathable film coverage

Breathable film material from front to back, top to bottom

The next generation of Prevention Plus is improved in every way. The gowns are exceptionally strong, yet comfortably light. They offer breathability and are film-reinforced throughout. Plus they're cut wider through the neck, chest and shoulders, for a roomy, comfortable fit. Our Prevention Plus fabric gives you a surgical gown that's technically advanced from end to end.

Breathable film reinforcement

Prevention Plus ensures you have film reinforcement all throughout your gown.

Cool comfort

Moisture vapor is smaller than the smallest known viruses. Using advanced technology, Prevention Plus allows moisture vapor to escape while preventing transmission of fluids.

Roomier fit

Prevention Plus gowns are cut wider through the neck, chest and shoulders for greater freedom of movement. They're available in extra-long and standard lengths to accommodate a variety of O.R. needs.



Item Number	Size	Non-Sterile Component
DYNJP2301P: 24/cs	Large (45", 114 cm)	19310PTM
DYNJP2302P: 24/cs	X-Large, Extra-Long (56", 142 cm)	21394PTM
DYNJP2303P: 22/cs	XX-Large, Extra-Long (56", 142 cm)	SPT-2303PTM
DYNJP2306P: 24/cs	Large (43", 109 cm)	SPT-2306PTM
DYNJP2307P: 24/cs	X-Large (47", 119 cm)	SPT-2307PTM
DYNJP2308P: 22/cs	XX-Large (49", 124 cm)	SPT-2308PTM
DYNJP2309P: 20/cs	XXXX-Large (49", 124 cm)	93043PTM



ADVANCING THE HEALTH OF HEALTHCARE™.

Medline exists to keep healthcare healthy. We're in an enviable position to do it. Privately held and 100% debt free, we have the freedom and flexibility to develop advanced medical and surgical products, and other innovative clinical and product solutions for your unique needs. Plus, the power and capital to deliver on our promises.

- » World-class network of more than 40 distribution centers
- » State-of-the-art manufacturing facilities—FDA QSR-compliant; ISO 13485-certified
- » Uncompromising quality – 550 dedicated quality professionals, inspection team of 100+ scientists and engineers at 17 quality locations



Know you're protected. Visit Medline.com for more information or contact your Medline representative today about a free trial!



Medline Industries, Inc. One Medline Place, Mundelein, IL 60060

Medline United States
1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-396-6996
medline.ca | canada@medline.com

Medline México
01-800-831-0898
medlinemexico.com | mexico@medline.com

FOLLOW US [f](#) [in](#) [t](#) [BLOG](#)