

# Advantage 200 5-Zone Mattress



As the trend toward home care continues to grow, so does the demand for products to help prevent and treat pressure injuries at home. The Advantage 200 5-Zone foam mattress provides a comfortable and cost-effective support surface to aid in the prevention and management of pressure injuries.

## Features

- Five distinct therapeutic zones for effective pressure redistribution
- 1-year warranty
- 300-lb. weight capacity
- Group 1 foam mattress
- Compression-packed for easy carrying and storage

Item No.	Description	HCPCS
MSCADV0280	Advantage 200 5-Zone Mattress 36 x 80 x 6" (91.4 x 203.2 x 15.24 cm)	E0184

Easy-to-clean nylon cover

Five distinct therapeutic zones for effective pressure redistribution

Non-skid bottom



Channel cuts create the 5 therapeutic zones of the mattress.



**Ease the pressure.** Contact your Medline Representative to request a sample or to place an order.