



TRIUMPH[®] NATURAL with ALOE

STERILE AND POWDER-FREE LATEX SURGICAL GLOVES

Many surgeons and operating room staff battle dry, irritated skin from constantly scrubbing with harsh surgical scrubs. In addition, the powder used in many surgical gloves as a donning agent can dry out your hands. In 1999, Medline changed the industry by patenting a process to bring the moisturizing properties of aloe vera to medical gloves.

Triumph Natural with Aloe Vera are latex, powder-free surgical gloves that will soothe and moisturize skin. They are designed for general surgical procedures where surgeons and staff that may have irritated skin.

Features and Benefits

- Ideal for general surgery and staff with dry, irritated skin
- Contains pure freeze-dried aloe vera to soothe and moisturize skin
- The lightly textured grip provides tactile sensitivity
- Powder-free and coated with a synthetic polymer coating to assist with damp hand donning
- ARC (anti-roll down cuff) was designed to reduce cuff roll-down common with most powder-free gloves
- Anatomical, curved former designed to enhance comfort and reduce hand fatigue

Specifications

- Natural Rubber Latex
- Synthetic Polymer Coating
- White
- Lightly Textured
- Standard Thickness 9.4 mil

Ordering Information

Item Number	Size
MSG2455	5.5
MSG2460	6.0
MSG2465	6.5
MSG2470	7.0
MSG2475	7.5
MSG2480	8.0
MSG2485	8.5
MSG2490	9.0

50 Pairs/Box, 200 Pairs/Case



TECHNICAL SHEET

TRIUMPH NATURAL with ALOE

Primary Material	Low-Protein Powder-Free Natural Rubber Latex with Synthetic Polymer Coating ≤2.0 milligrams/glove of powder in accordance with ASTM D6124 and ISO 21171 Caution: This product contains natural rubber latex which may cause allergic reactions, including anaphylactic responses.	
Donning Agent	Synthetic Polymer Coating (Inner surface coated for dry and damp hand donning) 100% Pure Freeze-Dried Aloe Vera Gel Coating (Moisturizes dry skin)	
Color	White (Provides contrast when using a dark colored underglove)	
Grip	Lightly Textured (Ideal for procedures requiring exceptional grip)	
Former (Mold) Design	Anatomical, curved former designed to enhance comfort and reduce hand fatigue	
Cuff Design	Tapered, beaded cuff design to prevent rolldown. Reinforced material prevents tearing.	
Chemical Additives (Accelerators)	Zinc Diethyldithiocarbamate (ZDEC) and Zinc Dibutyldithiocarbamate (ZDBC) Residual chemical levels below detectable level according to UPB/P/003a test method	
Leachable Protein (per EN455-3 using ASTM D5712) (Modified Lowry Protein Method)	Low Protein: <50 micrograms/gram of total extractable protein Caution: Safe use of these gloves by or on latex-sensitized individuals has not been established.	
Thickness (per ASTM D3577 ≥ 0.10 mm)	Finger Tip	0.24 mm
	Cuff	0.15 mm
Cuff Length (per EN455-2 ≥ 270 mm size 7.5)	300 mm	
Force @ Break Before Accelerated Aging (per EN455-2 ≥ 9 N)	Meets/Exceeds	
Force @ Break After Accelerated Aging (per EN455-2 ≥ 9 N, 7 days 70°C in an oven)	Meets/Exceeds	
Elongation @ Break Before Accelerated Aging (per ASTM D3577 ≥ 750%)	Meets/Exceeds	
Elongation @ Break After Accelerated Aging (per ASTM D3577 ≥ 560%, 7 days 70°C in an oven)	Meets/Exceeds	
Freedom from Holes (per EN 455 AQL 1.5)	0.65 AQL Final Inspection	
Viral Penetration	Tested and passed, in accordance with ASTM F 1671	
Chemical Resistance	The resistance to some chemicals has been assessed in accordance with EN 374-3 Results and recommendations for use with chemicals can be obtained on request	
Sterilization	Gamma Radiation, Sterility Assurance Level 10-6	
Expiration Date	35 Months from Date of Manufacture Manufacture and Expiration Dates are printed on packaging (YYYY-MM format)	
Packaging	Polyethylene peel pouch material protects product during transport and storage from moisture and ozone and prevents tearing when opening to maintain a sterile environment. Long pack configuration for easier aseptic donning.	
Storage Recommendations	Protect from freezing. Avoid excessive heat. Keep dry. Product should be shielded from direct sunlight, fluorescent lighting, X-rays, moisture and ozone. Do not store in temperatures above 40°C.	
Country of Origin	Malaysia	



Medline Industries, Inc.
One Medline Place
Mundelein, IL 60060

Medline United States
1-800-MEDLINE (633-5463)

www.medline.com | info@medline.com

Medline Canada
1-800-396-6996

www.medline.ca | canada@medline.com

Medline México
01-800-831-0898

www.medlinemexico.com | mexico@medline.com

Some products may not be available for sale in all countries. Please contact your local Medline sales representative, customer service or authorized Medline distributor for availability. We reserve the right to correct any errors that may occur within this sell sheet. ©2013 Medline Industries, Inc. Medline and Triumph are registered trademarks of Medline Industries, Inc., One Medline Place, Mundelein, IL.