



**EMPOWER**

**BEST PRACTICE.**

**Closed Suction Catheters**  
Better care. Optimal outcomes.

Innovating to inspire.

# Evidence-based care made easier.

You know your patients best, and your team works hard to ensure best practices are being followed to optimize outcomes. For closed suction, caregivers have been limited to choosing between products labeled for specific usage periods, which may conflict with following evidence-based guidelines regarding catheter change out policies.<sup>1</sup>

## Until now.

Introducing **ClearPro™** from Medline, your trusted partner in suction and tracheostomy care. ClearPro provides your facility the freedom to set—and follow—best practice protocols that can help optimize outcomes.

With an advanced cleaning chamber and enhanced sleeve, ClearPro empowers a better care experience for patients, while the Accessories Kit provides a better care experience for you.

**ClearPro: the power to set—and follow—best practice.**




1. AARC Evidence-Based Clinical Practice Guidelines. Care of the Ventilator Circuit and Its Relation to Ventilator-Associated Pneumonia. Respiratory Care. 2003;48(9):869-879. Available at: <https://www.aarc.org/wp-content/uploads/2014/08/09.03.0869.pdf>.

# New ClearPro™ Closed Suction Catheter System

Featuring a more advanced cleaning chamber and stronger sheath for better performance vs. competitors.

- » Labeled for facilities to determine change out policy
- » Does not contain DEHP
- » Isolated cleaning chamber for rapid catheter irrigation
- » Pressure wiper seal to help maintain PEEP while catheter is in use
- » Available in French sizes 10-16 with catheter depth markings and color-coded by ISO standards
- » Kitted with all necessary adapters and connectors to simplify SKU selection
- » High quality outer sleeve for a durable and tactile feel



# Medline Closed Suction Catheter Systems

## Adult

- » Labeled for facilities to determine change out policy
- » On/off valve allows suction catheter tip to be isolated and creates an enclosed cleaning chamber
- » Soft tip helps reduce damage to mucous membranes
- » Transparent connector for easy observation
- » High-quality, translucent EVA sleeve helps reduce cross-contamination and facilitates operation when advancing the catheter
- » One-way needleless irrigation
- » Catheter depth markings on all sizes
- » Color-coded locking valves

## Colored Thumb Valves



12 Fr



14 Fr



Also available with T-piece adapter



## Pediatric

- » Available with an elbow or Y-style modular patient adapter
- » Three adapter sizes to ensure compatibility
- » Labeled for facilities to determine change out policy
- » Compatible with all modes of ventilation
- » On/off valve allows suction catheter tip to be isolated and creates an enclosed cleaning chamber
- » Soft tip helps reduce damage to mucous membranes
- » Transparent connector for easy observation
- » High-quality, translucent EVA sleeve helps reduce cross-contamination and facilitates operation when advancing the catheter
- » One-way needleless irrigation
- » Catheter depth markings on all sizes
- » Color-coded locking valves

## Colored Thumb Valves



5 Fr



6 Fr



7 Fr



8 Fr

## Pediatric Adapters



Elbow adapters



Y-style adapters



# Ordering Information

## ClearPro Closed Suction Catheters

### Endotracheal

Item Number	Description	Pkg.	Thumb Valve Color
DYNCPDS10	ClearPro, Double Swivel Elbow, 10 French	20/cs	Black
DYNCPDS12	ClearPro, Double Swivel Elbow, 12 French	20/cs	White
DYNCPDS14	ClearPro, Double Swivel Elbow, 14 French	20/cs	Green
DYNCPDS16	ClearPro, Double Swivel Elbow, 16 French	20/cs	Orange
DYNCPTP10	ClearPro, T-Piece Adapter, 10 French	20/cs	Black
DYNCPTP12	ClearPro, T-Piece Adapter, 12 French	20/cs	White
DYNCPTP14	ClearPro, T-Piece Adapter, 14 French	20/cs	Green
DYNCPTP16	ClearPro, T-Piece Adapter, 16 French	20/cs	Orange

### Tracheostomy

Item Number	Description	Pkg.	Thumb Valve Color
DYNCPDS10T	ClearPro, Double Swivel Elbow, 10 French	20/cs	Black
DYNCPDS12T	ClearPro, Double Swivel Elbow, 12 French	20/cs	White
DYNCPDS14T	ClearPro, Double Swivel Elbow, 14 French	20/cs	Green
DYNCPDS16T	ClearPro, Double Swivel Elbow, 16 French	20/cs	Orange
DYNCPTP10T	ClearPro, T-Piece Adapter, 10 French	20/cs	Black
DYNCPTP12T	ClearPro, T-Piece Adapter, 12 French	20/cs	White
DYNCPTP14T	ClearPro, T-Piece Adapter, 14 French	20/cs	Green
DYNCPTP16T	ClearPro, T-Piece Adapter, 16 French	20/cs	Orange

## Medline Closed Suction Catheters

### Endotracheal

Item Number	Description	Pkg.	Thumb Valve Color
DYNCSDS12	Adult, Double Swivel Elbow, 12 French	10/CS	White
DYNCSDS14	Adult, Double Swivel Elbow, 14 French	10/CS	Green
DYNCSTP12	Adult, T-Piece Adapter, 12 French	10/CS	White
DYNCSTP14	Adult, T-Piece Adapter, 14 French	10/CS	Green

### Tracheostomy

Item Number	Description	Pkg.	Thumb Valve Color
DYNCSDS12T	Adult, Double Swivel Elbow, 12 French	10/cs	White
DYNCSDS14T	Adult, Double Swivel Elbow, 14 French	10/cs	Green
DYNCSTP12T	Adult, T-Piece Adapter, 12 French	10/cs	White
DYNCSTP14T	Adult, T-Piece Adapter, 14 French	10/cs	Green





## Medline Pediatric Closed Suction Catheters

Item Number	Description	Pkg.	Thumb Valve Color
DYNCPPE5	Pediatric CSC, Elbow Adapter, 5 French	10/cs	Black
DYNCPPE6	Pediatric CSC, Elbow Adapter, 6 French	10/cs	Light Green
DYNCPPE7	Pediatric CSC, Elbow Adapter, 7 French	10/cs	Purple
DYNCPPE8	Pediatric CSC, Elbow Adapter, 8 French	10/cs	Light Blue
DYNCPPY5	Pediatric CSC, Y-Connector, 5 French	10/cs	Black
DYNCPPY6	Pediatric CSC, Y-Connector, 6 French	10/cs	Light Green
DYNCPPY7	Pediatric CSC, Y-Connector, 7 French	10/cs	Purple
DYNCPPY8	Pediatric CSC, Y-Connector, 8 French	10/cs	Light Blue

**Empower optimal outcomes.**  
Contact your Medline Representative or call 1-800-MEDLINE to learn more about our complete offering of closed suction catheters and accessories.



Medline Industries, Inc. One Medline Place, Mundelein, IL 60060

FOLLOW US    

Medline United States  
1-800-MEDLINE (633-5463)  
medline.com | info@medline.com

Medline Canada  
1-800-396-6996  
medline.ca | canada@medline.com

Medline México  
01-800-831-0898  
medlinemexico.com | mexico@medline.com