

Tube Securement Devices and Dressings



Comfortable G-Tube securement.

With their skin-friendly hydrocolloid materials and long-lasting staying power, Medline's Tube Securement Devices and Dressings are a smart addition to your enteral feeding bundles.



Tube Securement Device

- A secure hold for feeding tubes, size 12-47 Fr
- Skin-friendly hydrocolloid adhesive to prevent injury; can be worn up to 7 days
- Soft and durable straps with no-slip grips for repeated positioning

Hydrocolloid Dressing

- Premium hydrocolloid adhesive provides a strong hold, yet is gentle enough to help reduce risks of skin tears and pressure injuries
- Low-profile for patient comfort
- Can be worn up to 7 days

Ordering Information

Item No.	Description	Pkg.
ENT51250NS	Tube Securement Device, With Strap, French Size 12-47	30/bx

Ordering Information

Item No.	Description	Pkg.
ENT54002NS	Hydrocolloid Dressing, Small	30/bx
ENT53002NS	Hydrocolloid Dressing, Medium	30/bx

Add tube securement devices and dressings to your enteral feeding bundles. To learn more, contact your Medline Representative.